

Product datasheet for RN207823

Spink6 (NM_001008874) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Spink6 (NM_001008874) Rat Untagged Clone
Tag: Tag Free
Symbol: Spink6
Synonyms: Spink511
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN207823 representing NM_001008874
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAAAGTAGCAGGTGTCTTTCTGCTCCTCTCTGGCTCTTCTCTGCTTTTTCTCAGGTGAGTTAGTC
 AAGGAGGACAGGACAAAAGAGGCTGGATACAAGATCAGAGGGCCGATTCCAGGAAAGGTGTCTTAAG
 GCATCGTCTCTTTCAAATTAATTGTGGTGAGTTCCGTGATCCCAAGGTCTTCTGCACTCGGGAATCTGAC
 CCACTCTGTGGCTCTGATGGCCAGACATATGGGAATAAATGTGCCTTCTGCAAGGCATTGGAGAAAAGTT
 CTGGAAAGATCAACCTGAAGCATCGTGGAAATGCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul

ACCN: NM_001008874

Insert Size: 318 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001008874.1](#), [NP_001008874.1](#)

RefSeq Size: 318 bp

RefSeq ORF: 318 bp

Locus ID: 408235

UniProt ID: [Q6IE47](#)

Cytogenetics: Xq21

Gene Summary: predicted to act as a serine protease inhibitor [RGD, Feb 2006]