

Product datasheet for **RN207846**

Fabp6 (NM_017098) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fabp6 (NM_017098) Rat Untagged Clone
Tag: Tag Free
Symbol: Fabp6
Synonyms: IBABP
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN207846 representing NM_017098
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCCTTACC GGCAAATACGAATTTGAAAGTGAGAAGAATTACGACGAGTTCATGAAGCGCCTGGGTC
 TTCCAGAAGATGTGATTGAAAGGGGACGTAACCTCAAGATCATCACAGAGGTCCAGCAGGATGGAGAGAA
 CTTCCACTGGTCCAGTCTTACTCTGGGGCAACATCATGAGCAACAAGTTCACCATTGGCAAAGAATGT
 GAAATGCAGACCATGGGGGGCAAGAAGTTCAAGGCAACCGTGAAGATGGAGGGTGGCAAGGTGGTGGCAG
 ACTTCCCAACTATCACCAGACTTCGGAGGTCTGGGTGACAAGTTGGTGGAGATCTCCACCATCGGGGA
 TGTGACCTACGAGCGTGAAGCAAGAGGGTGGCT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_017098
Insert Size: 387 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_017098.1](#), [NP_058794.1](#)

RefSeq Size: 498 bp

RefSeq ORF: 387 bp

Locus ID: 25440

UniProt ID: [P80020](#)

Cytogenetics: 10q21

Gene Summary: may mediate metabolism of steroids in steroid hormone-producing tissues [RGD, Feb 2006]