

Product datasheet for **RN207897**

Mfsd13a (NM_001014030) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Mfsd13a (NM_001014030) Rat Untagged Clone
Tag: Tag Free
Symbol: Mfsd13a
Synonyms: RGD1309313; Tmem180
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN207897 representing NM_001014030
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGC**C

ATGGGGCTAGACTGGCCCCAGCCTGGCTGTTGGGCTGCCACAGCCGTGGTGTATGGCTCCTTGCTT
TCTTCACCACCATTCTGCATAATGTGTTTCTGTGCTACTACGTGGACACCTTTGTCTCGGTGATAAAGT
CAACAAAGTGTCTTCTGGGTTGGAGAGACTGTATTTCTCTCTGGAATAGCTTCAATGACCCCTCTTC
GGTTGGCTCAGTGACCGCGGTTCTCAGCTCCCAGCCCGGTCAGGAGCAGGGCTCTCCTCCAGAGATG
TGGTGTGACCCGGGTGCGGGCTCTGGGTTGGCAGGGCCACTGCTGGCACTGTCCTTCTGGCATTCTG
GGTGCCCTGGGCCCTGCTGGCTTACAGTTCTTACTGTGCCTGTGCCTTATGATGGCTTCTGACACTC
GTGGACCTCCACCATCATGCCTTGTGGCCGACCTGGCGCTCTCGGCCATGACCGCACCCACCTCAACT
TTTACTGCTCCCTCTTCAAGTGCAGCTGGCGCCATGTCTGTCTTTGCCTCTACGCCTTTTGAACAAGGA
GGACTTCTCTCTTCCGTGCCTTCTGTGTGGTACTAGCTGTGGCTCTGGGCTGGGCTTTTGGGGCC
ACACAATTGTTGAAGCGGCGATTGAGGCGACCAGAAAGGACAGGGGTTGCCGGCCTGGTCATGGATG
GTGGCGTGTGTAAGAGGAGCCGCTGGTGGTGGTGGGAGAGCAGCAGCATTACCTTGGGTGAGTATCT
CCGGCAACTGGCACACCACAAAACCTTCTGTGGTTCGTGGGCATGGACCTGGTACAGGTGTTTCACTGC
CACTTCAACAGTAACTTCTTCCCTCTTCTGGAGCACTTGTGTGGACACATCTCTCTCCACTG
GCTCCTTCTGTTGGCAACCGTGTCTTACCGAGGGTACCTGTAAGCTG**CTGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001014030
Insert Size: 963 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001014030.1</u> , <u>NP_001014052.1</u>
RefSeq Size:	2064 bp
RefSeq ORF:	963 bp
Locus ID:	309454
Cytogenetics:	1q54