

Product datasheet for RN207912

Lhx2 (NM_001106571) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lhx2 (NM_001106571) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Lhx2
Synonyms:	Lh-2; LH2A
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN207912 representing NM_001106571 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGTTCCACAGTCTGTGGGCCCCGAGGTGCACGGGGTCATCGACGAGATGGACCGCAGGGCCAAGA
GCGAGGCTCCGGCCATCAGCTCCGCCATCGACCGCGGACACGGAGACGACCATGCCGTCATCAGCAG
TGACCGGGCAGCGTTGTGTGCGGGCTGTGGGGCAAGATCTCTGACCGCTACTATTTGCTGCCAGTACAG
AAACAATGGCACATGCGTTGCCTCAAGTGTGAATGCAAGCTCAACCTGGAGTCGGAACCTCACCTGCT
TCAGCAAGGATGGCAGCATCTACTGCAAAGAAGACTACTACAGGCGGTTCTCTGTGCAGCGCTGCGCCCG
CTGCCACCTGGGCATCTCGGCCTCAGAGATGGTGATGCGTGCTCGGGACTTGTTTATCACCTCAACTGC
TTCACATGCACCACGTGTAACAAGATGCTGACGACCGGCGACCACTTCGGCATGAAGGACAGCCTGGTCT
ATTGCCGTTGCACTTCGAGGCACTGCTGCAAGGCGAATACCCGGCACACTTTAACCATGCCGACGTGGC
AGCGGGCGCAGCGGCAGCCGAGCAGCTAAGAGTGCAGGATTGGGCTCAGCTGGGCTAACCCGCTGGGT
CTTCCCTACTACAACGGCGTGGCACTGTGCAGAAGGGGAGGCCGAGAAAGCGCAAAAGTCCAGGACCCG
GGCAGATCTGGCAGCCTACAACGCCGCTGAGCTGTAACGAGAACGACGCGGAACACCTGGATCGCA
CCAGCCCTACCCTAGCAGCCAAAAGACAAAGCGCATGCGCACCTCCTTCAAGCACCACGCTTCGGACA
ATGAAGTCTTACTTTGCCATTAACCACAATCCCGATGCCAAGGACTTGAAGCAGCTCGCGCAAAAGACCG
GCCTCACCAAGAGAGTCTCCAGGTCTGGTTTTCAGAATGCCCGGGCAAGTTCAGGCGCAACCTTTTACG
GCAGGAAAACACTGGTGTGGACAAGACTTCGGATGCCACGCTGCAGACAGGGACACCGTCAGGGCCCGCC
TCGGAGCTGTCCAACGCCTCGCTCAGCCCTCCAGCACTCTACAACCCTCACAGACTTGACTAGCCCCA
CCCTGCCGACTGTGACGTAGTCTTAACCTCTGTGCTGGCAACCTGGAGGGCCATGAGCCCCACAGCCC
TTCACAAACGACTCTTACCAACCTTTTCTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001106571
Insert Size:	1221 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001106571.1, NP_001100041.1</u>
RefSeq Size:	1903 bp
RefSeq ORF:	1221 bp
Locus ID:	296706
Cytogenetics:	3q12
Gene Summary:	LIM domain transcription factor; involved in lymphoid and neural differentiation [RGD, Feb 2006]