

Product datasheet for **RN208022**

Foxs1 (NM_001012091) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Foxs1 (NM_001012091) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Foxs1
Synonyms:	Fkhl18
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN208022 representing NM_001012091 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCAGCAGCCACCCACAGCAGAATCCTTGCTCCTTCTACAGAGCCCAGCAAGCCCCCTTACAGCTACA
TAGCCCTGATTGCCATGGCTATCCAGAGTTCACCGGGTCAGCGCGCCACCCTCAGTGGCATCTACCGCTA
CATCATGGGCCGCTTTGCATTTTACCGGCACAACCGGCCAGGCTGGCAAAACAGCATCCGCCATAACCTG
TCTCTCAACGAGTGCTTCGTCAAGGTGCCTCGTGATGACCGCAAGCCAGGCAAGGGCAGCTACTGGACCC
TGGATCCAGACTGCCACGACATGTTCCAGCACGGCAGCTTCTCCGTCGCCGAGACGCTTACCAAGCG
CACAGGTGCCAGGGACCAAGGGCCCTGTCAAGGCAGACCACAGACCCCTTAGAGCCACCAGCCCGGAT
CAAGGAGCCCCAACACCACAACCTGGCAGACTGTGCCATTCCCACCAGAGGTGCCAAACCCCAAGGGCT
TTGGGGGTCTGATGGGGTCATTGCCAGCCAACATGTGCCCAACAACCAAGTACACCCAGGCCCCAGCTGCC
CACTGGACCCAAAGACATGTGCTCGGCAAAGTCTGGAGGTCCAAGGGAGCTTTCTGAAGCCACCTCGCCC
TCACCATGCCCGGCAATTCGGCTTCTCATCTGCCTTCTCTGAAGCTGAGAGTCTCGGCAAGGCCCCACAC
CCAGTGTGGTCCGGAGAGCATCGGGAGCAGCTACCAGTGCCGGATGCAGACGCTGAATTTCTGTATGGG
GGCTGACCCAGGCCTTGAACACCTTTGACCTCTGCAGTCGCCACTCCTGGATCCTCAACACCTTCGGCC
TCTCACCAGGCTCCACTGCCCTGCCAGCTGACTCCAAGGAACCGTGGGTGCCTGGGGCTTCCCTGTCC
AGGGAGGTTCTGGCTATCCACTGGGATTACCCCTTGTCTTTACCGGACACCGGGAATGTTCTTCTCGA
GTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul



[View online »](#)

ACCN:	NM_001012091
Insert Size:	984 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001012091.1</u> , <u>NP_001012091.1</u>
RefSeq Size:	1300 bp
RefSeq ORF:	984 bp
Locus ID:	311547
Cytogenetics:	3q41