

## Product datasheet for **RN208046**

### Dusp11 (NM\_001025650) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Dusp11 (NM\_001025650) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Dusp11  
**Synonyms:** MGC112701  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN208046 representing NM\_001025650  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGAACCACTGGCATTATGGCCGATACAGTCGGGGCCGGGACTTTACTGCCCGTGCCCCACCCAAGAAGA  
AGGGCAAAAACCAGATCCCGAAAGGTGAAAGACTATCTCCAGTTGGGCAACGGATGCCTGGGACTCG  
TTTCATTGCTTTCAAAGTTCCTTTGCAAAAAAATTTGAAGCAAAGCTTATGCCAGAAGAATGCTTTTCT  
CCTTTGGATCTCTTAATAAAATCAAGAACAAAATGAAGAACTTGGATTGATAATCGATTAACATATA  
CTCAGCGCTATTATAAAGTAGAGGATTTGCCAAAAACTATTTCTTACATAAAAAATCTTACAGTTGGTCA  
TCAGGTACCTGACAGTGGCACTATTTTTAGTTCAAAGTGCAGTTAAAGAGTTTTTAAAGAGAAAATAAA  
ACAATGACAACTTATTGGTGTCCACTGTACCCATGGCTTAAATAGGACTGGCTACCTCATTGCGAGAT  
ATTTGATTGATGTAGAAGGCATGAGGCCAGATGATGCAATTGAATTGTTCAACAGGTGCCGGGACATTG  
CATAGAAAGGCAGAACTACATCGAAAACCTGCAGAAACGTCGTGTGAGAAAGAACAGAAATGCTAGTGCA  
TCTAGATCAGGTGGTCTTGAAGACTCAGCACATTTACAGAACAAGTCCACACCACTAACAAACAGTGGA  
ACAAAGGACCTAAGAAAAGCCGACGTGGTGGTCACTGGAATCGTCCCAGCATGTTCCAGCCAGCGCTC  
AGCGTACTCATTAGAAAATGGTCACAAAATCAGAGTGTACCAGAGAGGCTTCGTCCCTCCTCCCGGT  
CCAGCTGGAGAGGACTATTCACAGAGAAGATTCTTCTGGAGCATGAGACCAAATGGCAGTCAAGCAACCC  
ATCATAAGAAATGGATTGCTGCTAGTTACCAGAGACCTTTCTACCCAGCCTCCTGGGAATGGAATGTGTG  
A

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001025650



[View online »](#)

<b>Insert Size:</b>	981 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001025650.1</a></u> , <u><a href="#">NP_001020821.1</a></u>
<b>RefSeq Size:</b>	1155 bp
<b>RefSeq ORF:</b>	981 bp
<b>Locus ID:</b>	297412
<b>UniProt ID:</b>	<u><a href="#">Q4KM79</a></u>
<b>Cytogenetics:</b>	4q34
<b>Gene Summary:</b>	Possesses RNA 5'-triphosphatase and diphosphatase activities, but displays a poor protein-tyrosine phosphatase activity. In addition, has phosphatase activity with ATP, ADP and O-methylfluorescein phosphate (in vitro). Binds to RNA. May participate in nuclear mRNA metabolism.[UniProtKB/Swiss-Prot Function]