

## Product datasheet for **RN208048**

### Mtx2 (NM\_001008286) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Mtx2 (NM\_001008286) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Mtx2  
**Synonyms:** MGC106000  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN208048 representing NM\_001008286  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGTCTCTGGTGGCGGAAGCGTTCGTCTCCAGATCGCAGCTACAGAACCATGGCCTGAAAATGCAACCT  
TATATCAGCAACTGAGAGGGGAGCAAATTTTGCTTTCTGACAATGCAGCATCTCTTGCTGTGCAGACCTT  
TTTGCAGATGTGAATCTGCCTGTCAAAGTGGTATGCAGGGCAAATGCAGAATATATGTCTCCATCTGGA  
AAAGTACCTTTTATTCATGTGGGAAATCAAGTAGTGCAGAACTGGCCCAATAATCCAATTTGTTAAAG  
CCAAGGGCCATTCTTCTAGTGATGGGTTGGATGAAGTCCAAAAAGCAGAAATGAAGGCTTACATGGAATT  
AGTCAACAACATGCTGCTGACTGCAGAGTTGTATCTCCAGTGGTGTGATGAAGCTACAGTAGGGGAGATC  
ACTCTTGCTAGGTATGGGTCTCCTTATCCTTGGCCTCTGAACCTCATTCTGACCTACCAGAAGCAGTGTG  
AAGTCAAGCGTAAGATGAAAGCCATTGGATGGGGTAACAAGACCCTGGACCAGGTCTTAGAAGACGTAGA  
CCGGTGTGTCAAGCCCTTCCCAAAGACTGGGAACACAACCTTACTTCTCGATAAGCAACCTACTGAA  
CTTGATGCTCTGGTTTTTGGCCATTTGTACACCATTCTTACCACACAGTTGACCAGTGTGAACCTTTGTG  
AGAAGGTGAAAACTATAGCAACCTCCTTGCTTCTGTAGAAGAATTGAACAGGACTATTTTGAAGACCG  
GGTAAAGGCAGGTTGTCT**TAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001008286  
**Insert Size:** 792 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001008286.1</a></u> , <u><a href="#">NP_001008287.1</a></u>
<b>RefSeq Size:</b>	1235 bp
<b>RefSeq ORF:</b>	792 bp
<b>Locus ID:</b>	288150
<b>Cytogenetics:</b>	3q23