

Product datasheet for **RN208083**

Glo1 (NM_207594) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Glo1 (NM_207594) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Glo1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN208083 representing NM_207594 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGCAGAGCCACAGCCAGCGTCCAGTGGTCTTACTGATGAGGCCGCCCTCAGTTGCTGCTCTGACCCGG
ATCCAAGCACCAAGGATTTTCTACTGCAGCAAACGATGCTGAGAATTAAGGACCCCTAAGAAGTCCCTGGA
TTTTTATACGAGGGTTCTTGACTGACGCTTCTCCAGAAGCTTGACTTCCCTTCTATGAAGTCTCGCTC
TATTTCTTAGCGTACGAGGATAAGAACGATATCCCAAGGACAAGACGGAGAGGACAGCGTGGGCTTTTT
CCAGAAAGGCTACCCTGGAGCTCACACAACTGGGGCACTGAAGATGACGAGACGCAGAGTTACCACAA
CGGCAACTCGGACCCTCGGGGATTTGGCCACATTGGGATTGCTGTTCTGATGTCTACGAGGCTGTAAG
AGATTTGAAGAACTAGGTGTCAAGTTTGTGAAGAAGCCTGATGATGGGAAAATGAAAGGTCTGGCGTTCC
TTCAAGATCCTGATGGCTACTGGATTGAAATTCTGAATCCTAACAAAATGGCAACAATTATT**AG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_207594
Insert Size:	555 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_207594.2](#), [NP_997477.1](#)

RefSeq Size: 1079 bp

RefSeq ORF: 555 bp

Locus ID: 294320

UniProt ID: [Q6P7Q4](#)

Cytogenetics: 20p12

Gene Summary: predicted to have lactoylglutathione lyase activity [RGD, Feb 2006]