

## Product datasheet for **RN208185**

### Adprhl2 (NM\_001108680) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Adprhl2 (NM\_001108680) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Adprhl2  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >RN208185 representing NM\_001108680  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCGGCGGGCGGGCTGCGGCAGCGGCAGCGATGACGGCGGGCGGGCGGGGAGCTGGTGCGGCC  
 GCTCCCTCTCGCGCTCCGAGGCTGCCTGGCGGGCGCGCTGCTGGGAGACTGCGTGGGCGCTGTCTTCGA  
 GGCACACGATACCGTCAGCCTGACGTCTGTCTGAGTACGTTTCAGAGCCTGGAGCCAGACCCGGGCAGC  
 CCGGGCAGCGCGGGACAGAAACATTGTACTACACAGATGACACCGCCATGACCAGGGCCATGGTGCAGT  
 CCCTGCTGGCCAAGGAGGCTTTTCGACGAGGTGGACATGGCTCACAGTTTGCCCAAGAATAAAGAAGGA  
 CCTGACAGAGGGTATGGGGCCGAGTCATCACTGTCTTCAAGAACTCCTGAATCCAAAGTGCCGTGAT  
 GTGATGAACCTGCCCGGCCAGTTCAATGAAAGGGCTCCTATGGCAATGGGGAGCCATGCGGGTAG  
 CTGGCATCCCGCTGGCCTATCACAGTGTCCAAGATGTGCAGAAGTTTGCCCGCTCTCAGCTCAGCTGAC  
 CCATGCCTCTTCCCTGGTTACAACGGTGCCATCTTGCAGGCCCTGGCTGTGCACCTAGCGCTGCAGGGT  
 GTATCATCCAGTGAGCACTTCTCGAGCAGCTTCTGGCCACATGGAGGAGCTGGAAGGTGATGCCAGT  
 CAGTCTTGGATGCCAAGGAGTTGGGTATGGAGGAGCGCCCTACTCCAGCAGGCTGAAGAAGGTGGGAGA  
 GCTGCTGGACCAGGACCGGTGAGCCGAGAGGAAGTGGTGTCTGAGCTAGGAATGGAATCGCTGCCTTT  
 GAGTCTGTGCCACCGCCATCTACTGCTTCTGCGCTGCATGGAGCCTGATCCTGAGATCCCTCCACCT  
 TCAACAGTCTCCAGAGGACTCTCATCTACTCCATCTCTTGGTGGGACACAGACACCATAGCCACCAT  
 GGCTGGGGCCATTGCTGGAGCTTACTATGGGGTGAACAGGTGCCGGAGAGCTGGCAGCAGAGCTGTGAG  
 GGCTATGAGGAGACAGATGTCCTGGCCAGAGCCTACACCGAGTCTTCCAAGAGAGTCCCTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001108680



<b>Insert Size:</b>	1113 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001108680.1</a></u> , <u><a href="#">NP_001102150.1</a></u>
<b>RefSeq Size:</b>	1666 bp
<b>RefSeq ORF:</b>	1113 bp
<b>Locus ID:</b>	362600
<b>Cytogenetics:</b>	5q36