

Product datasheet for **RN208216**

Abhd8 (NM_001107301) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Abhd8 (NM_001107301) Rat Untagged Clone
Tag: Tag Free
Symbol: Abhd8
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN208216 representing NM_001107301
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACCTCAGGACACACCATGCTGACCGGCGTGACAGATGGCTTCTTTTGCTGCCTGCTGGGGACTCCCC
CTAACCCCGTGCGGCCGCTGGAGAGCGTCGAGTCCAGTGACGGCTACACGTTTCGTGGAGGTCAAACCAGG
CCGGGTGCTGCGGGTAAAGCATGCGTGCCCGCGCCCATCCCCTACCGCCACCCCTCCTCCTGAAGAT
AATCCAGGGGTCAAGACGGGCCTTGTGCGCTGCCAGCGCGAATCACAGTGTACCGAAATGGGAGGCTGG
TGGTGGAGAACCCTGGGCCGAGCGCCTCGTGCCGACCTGCAAGGCCGACGTGGCTCCGGGGACCCGCTGG
CTCACTTGAGGTAGAGTTGGCTGAGCCCGAGGAGGCGACACCCGGGCCACCCGGGCAGCGGGCGACGC
CGCCGTCCACGGCGTCCGAAGAGGACCATCCACATCGACTGCGAGCAGCGCATACCAGCTGCAAGGGCG
CGCAGGCTGATGTGGTGTCTTTCATCCATGGAGTGGTGGCTCACTCGCCATCTGGAAGGAGCAGCT
GGACTTCTTTGTGCGACTCGGTTACGAGGTGGTGGCACCTGACCTGGCCGGCCATGGAGCCAGCTCTGCG
CCACAGGTTGCAGCCGCCTACACCTTTACGCATTGGCAGAGGACATGCGGGCCATTTTCACGCGCTATG
CCAAGAAGCGCAACGTGCTCATTGGTCATTCTTACGGGGTCTCCTTCTGTACTTTCTGGCCATGAGTA
CCAGATCTTGTACACAAGGTGATGATCAACGGTGGAGGCCACCGCACTGGAGCCAGTCTCTGC
TCCATCTTCAACATGCCACGTGTGTCCTGCACTGCCTGTACCTTGCCTGGCCTGGAGTCTCTCAAGG
CCGTTTTGCTCGCAAGGGGCCAAGGAAAAGCAGCTGCTGAAGGAGGGCAATGCCTTCAATGTGTCCTC
TTCGTGCTGCGGGCCATGATGAGCGGCCAGTACTGGCCTGAAGGTGACGAGGTACCATGCGGAGCTG
ACAGTTCCCGTGCTTCTGGTGCACGGCATGCACGACAAGTTTGTGCCAGTGAAGAGGACCAACGCATGG
CTGAGATCCTGCTGCTGGCCTTCTCAAGCTCATCGAGGAAGGCAGCCACATGGTGTGCTGGAGTGTC
TGAGACAGTCAACACACTGCTACACGAGTTCCTGCTGTGGGAGCCAGAACCTGAGGCCGAGCCAAAGCTC
GAGCCAAAGTCCAAGCCACAGCCGCTGCAGCTGCAGTCAGACCCGGCTCCAGGGGAGGAGAAG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



ACCN:	NM_001107301
Insert Size:	1326 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001107301.1, NP_001100771.1</u>
RefSeq Size:	1982 bp
RefSeq ORF:	1326 bp
Locus ID:	306338
Cytogenetics:	16p14