

Product datasheet for **RN208252**

Rps3a (NM_017153) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rps3a (NM_017153) Rat Untagged Clone
Tag: Tag Free
Symbol: Rps3a
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN208252 representing NM_017153
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGTGCGGAAGAACAAGCGCCTGACGAAAGGTGGCAAGAAGGGAGCCAAGAAGAAAGTGGTTGATC
CATTTTCTAAGAAAGACTGGTATGATGTGAAAGCTCCAGCAATGTTCAATATTAGGAACATCGGGAAAAC
ACTAGTTACAAGGACTCAAGGAACAAAATTGCATCAGATGGCCTCAAGGTCGTGTGTTTGAAGTGAAT
CTTGACAGATCTACAGAACGATGAAGTTGCATTTAGAAAGTTCAAGCTAATTACTGAGGATGTTTCAAGGCA
AAAAGTGTCTGACTAACTTCCATGGCATGGATCTTACCCGGGACAAAATGTGCTCCATGGTCAAAAAGTG
GCAGACCATGATTGAAGCTCATGTTGATGTGAAGACGACCGATGGGTATTTGCTCCGACTTTTCTGTGTT
GGTTTCACTAAGAAACGTAACAACCAGATTAGAAAGACATCCTATGCGCAGCACCAGCAGGTCCGCCAGA
TCCGCAAGAAGATGATGGAGATCATGACCCGAGAAGTGCAGACAAATGACTTGAAGGAAGTAGTTAATAA
ACTGATCCAGACAGCATTGGGAAAGACATAGAGAAGGCTTGCCAGTCCATCTATCCTCTCCATGACGTT
TTTGTAGGAAAGTCAAATGCTGAAGAAACCAAGTTTGAATTGGGAAAACATGGAGCTCCATGGTG
AAGGCGGTAGTTCTGGCAAAACAACCTGGTATGAGACAGGTGCTAAGGTTGAACGAGCTGATGGCTATGA
ACCACCGTCCAAGAATCCGTT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul
ACCN: NM_017153
Insert Size: 795 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_017153.1 , NP_058849.1
RefSeq Size:	880 bp
RefSeq ORF:	795 bp
Locus ID:	29288
UniProt ID:	P49242
Cytogenetics:	2q34
Gene Summary:	40S ribosomal subunit protein [RGD, Feb 2006]