

Product datasheet for **RN208330**

Ccl3 (NM_013025) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ccl3 (NM_013025) Rat Untagged Clone
Tag: Tag Free
Symbol: Ccl3
Synonyms: MIP-1a; Scya3
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN208330 representing NM_013025
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAAGGTCTCCACCGCTGCCCTTGCTGTTCTTCTCTGCACCATGGCGCTCTGGAACGAAGTCTTCTCAG
CGCCATATGGAGCTGACACCCCGACTGCCTGCTGCTTCTCCTATGGACGGCAAATCCACGAAAATTCAT
TGCTGACTATTTGAGACCAGCAGCCTTTGCTCCAGCCGGGTGCATTTTCTGACCAAGAGAAACCGG
CAGATCTGCGCTGACCCAAAGAGACCTGGGTCCAAGAATACATCACTGAGCTGGAACAAATGCC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul
ACCN: NM_013025
Insert Size: 279 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_013025.2](#), [NP_037157.2](#)

RefSeq Size: 755 bp

RefSeq ORF: 279 bp

Locus ID: 25542

UniProt ID: [P50229](#)

Cytogenetics: 10q26

Gene Summary: mediates monocyte and neutrophil chemotaxis; may play a role in the pathogenesis of acute lung injury [RGD, Feb 2006]