

## Product datasheet for **RN208442**

### **Pdss2 (NM\_001014249) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Pdss2 (NM\_001014249) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Pdss2  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN208442 representing NM\_001014249  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAGCCTCCGGCAGCTGCTGTTGCGCTTGTCCGGTTATCTCGGGGCTTCTGGTCCCCAAATCGCCACT  
GGTGGTACTTCAGATCTCTCGACACCATCTCCTCGGTAGGCTCCTGGCGTGGTCCGCTCCAGGTCTCC  
AGCCCATTTGGAACAGGTGGTGTGACAGGCAGAGAAGATCGTGGGTTACCCCGCATCCTTCATGAGCCTC  
CGCTGCCTGCTGAGCGACGAGCTCAGCAACGTCGCTATGCAGGTGCGGAAGCTGGTGGGGACTCAGCACC  
CACTGCTTACCACTGCCAGGGGCTTCGTGCATGACAGCAGACATAACCTACAGCTGCGGGGCTGGTGGT  
GCTGCTCATTCCAAAGCGGGGCCAGCACTAGGAACTCTTCCAGTCAGAACTACGACATGGTCAGT  
GGGATATACTCATGTCAAAGAAATTTGGCAGAGATCACAGAACTTATACATACTGCTCTCCTGGTGCATC  
GTGGGATAGTAAACCTAAGCGAATTACAGTCTTCTGATGGACCACTGAAAGACATGAAGTTTGGAAACAA  
GATAGCTGTCCTGAGTGGAGACTTCTCCTAGCAAACGCCTGTAATGGACTAGCTCTTCTACAGAACACC  
AAGTTTGTGGAGCTCCTAGCAAGTCTGATGGACTTGGTCAAGGCATCTACCAGGAAAACCTGCTT  
CCACCCAGGGAAATCCTATCCAGATGACATCAGAATATCGACCTGGAAGGAGCAGACTTTCCTATCCCA  
CTGTGCCTTGTAGCGAAGAGCTGCAAGCCGCCATGGAGTTAGCGAAGCATGACGCTGCGGTCCAAGAC  
ATGGCGTTTCAATACGGAAGCACATGGCCATGAGCCACAAGATAAATCTGACCTCCAGCCTTTTATTA  
AAGACAAGGCCAGTGACTCTAAGACTTTTAACTAACTCGGCACCCGTAAGTCTTACATCAGGAGTTTCT  
CGGGAGAGATTTGTGGATTAAGCAGATTGGAGAGGCTCAAGAAAAAGGAAGTTGAACTACACTAAGTTG  
CGGAAACAATCAAAGCTGGCAAAGGTGTGACTTCAGCTATTGACCTGTGTGCTTACCATGGAAACAAGG  
CACTGGAGGCCCTGGAGAGCTTCCCTCCCTCAGAGGCCAGATCGGCTTTAGAAAAATTGTGTTTGGCGT  
GACGAGATTTTCA**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI



[View online »](#)

<b>ACCN:</b>	NM_001014249
<b>Insert Size:</b>	1206 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001014249.1</a></u> , <u><a href="#">NP_001014271.1</a></u>
<b>RefSeq Size:</b>	2100 bp
<b>RefSeq ORF:</b>	1206 bp
<b>Locus ID:</b>	365592
<b>UniProt ID:</b>	<u><a href="#">Q5U2R1</a></u>
<b>Cytogenetics:</b>	20q13
<b>Gene Summary:</b>	Supplies decaprenyl diphosphate, the precursor for the side chain of the isoprenoid quinones ubiquinone-10.[UniProtKB/Swiss-Prot Function]