

## Product datasheet for **RN208498**

### Spred1 (NM\_001047089) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Spred1 (NM\_001047089) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Spred1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN208498 representing NM\_001047089  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGAGCGAGGAGACGGCGACTTCTGACAACGATAATAGTTATGCACGAGTGCGAGCTGTGGTGATGACCA  
GAGATGACTCAAGTGGTGGATGGCTACCACTCGGAGGGAGCGGACTGAGCATTGCTACTGTCTTCAGAGT  
CCCTCAGCAGGAAGAGAATGGCTGTGACAGCTTTTATATCCGTGGAGAGAGACTCAGAGACAAAATGGTG  
GTTTTGGAATGTATGCTTAAGAAAGACCTCATCTATAATAAGGTCACTCCCACTTTTCACCACTGGAAGA  
TCGATGACAAGAAGTTTGGCCTTACCTTTCAGAGTCTGCTGATGCCAGGGCTTTTGATCGAGGCATTTCG  
AAGAGCTATAGAAGATATTTCTCAAGGGTGCCAGCATCAAAAAGTGAAGTGAAGGAGGAGATGATGAT  
TTACAAACAAGTGAAGGACACTTCCCGTTCCCTGGTAAAGATCACTTTTTTCAGCAAGAGACAGTTG  
TTACCTGTGATCCTTACAGAAGCTCAGACATAAGACCATTACCCTTTGAAGATCTGAATGCCAGAAGAGT  
CTACTTGCAAAGCCAAGTGAAGCAGATACCACTTTGTGAGCAAGGCTAGACATACAGAGCCGAAGTATG  
GAATATGTACAGCGGCAGATATCTAAGGAATGTGGGAGCCTCAAGTCCCAAGCTAGGGTCCCTTGAAGT  
CTATCAGACATGTCAGCTTTCAAGATGAGGATGAGATTGTTAGAATAAATCCTCGAGATATCTTAATACG  
TCGGTATGCAGACTACAGACATCCGGACATGTGGAAAAATGACTTGGAGAGAGATGATACTGACTCCAGT  
GTTCCATTTTCTAAACAAGACAGTAAAAAGTCAAGTATTTCTACCACTGTGGGGATGAGACTAAGTTAA  
GTTCTCTGAAAGACTCTGTGGTATTTAAGACACAGCCTCCCTCATTAAAAATTAAGTCAAAACGAAGAAA  
AGAGGATGGGGAACGTTCTCGTGCCTGTACTGCCAGGAAAGGTTAATCATGAAGAAAATGGCAGGGGG  
AAATGCCAGGATGCTCCAGATCCTGTTAAAAGATGCATATATCAAGTTAGTTGCATGCTCTGTGCCGAGA  
GCATGCTGATCATTGTATGTGACACTCAGAGGGAGACTTTTCTGACCCTTGTTCGTGTGACACTAGCGA  
TGACAAGTTCTGCCTAAGGTGGTTAGCCCTGGTAGCTCTGTCTTTCAATTGTACCGTGTATGTGCTGCTAC  
GTCCCTTTGAGAATGTCCATCGCTGTGGTGGGATGTGGTGTGGTGGGAAACACAAGGCTGCTG  
GGTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAAAGTCACTCAGAAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_001047089
<b>Insert Size:</b>	1335 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001047089.1</a>, <a href="#">NP_001040554.1</a></u>
<b>RefSeq Size:</b>	1335 bp
<b>RefSeq ORF:</b>	1335 bp
<b>Locus ID:</b>	296072
<b>Cytogenetics:</b>	3q35