

Product datasheet for **RN208561**

Ctif (NM_001108891) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ctif (NM_001108891) Rat Untagged Clone
Tag: Tag Free
Symbol: Ctif
Synonyms: Gm672
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN208561 representing NM_001108891
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAAGCATCCGCCAGGAGACAAGGGAGAAGCGGGCAGCCACCGCAACGCCAAAGAGACCGTGACCGTTG
AGAACCCCAAACCTGAAGACGGCCAGGGGACACCGGACACAGTGGCCTTGAGCCCCCTGCAGCCCTGA
CACTCTGACCCCAACAGCTTCCGAGAGGCCAACCCCTCAGCTGCCAGGGGTTTCCAGAGGCTGAGACAAAG
CAAAAAGACACTGCTTCCCGAGCGTCTTGGGGAGCGGCCAAAATCACCTGCTCCAGTCTTCCAAAG
ACAGGCTGAGGCGAAGGCTGAAAGAAAAGTACCGGACGAAGTAGCAGTAGAGGCCAGCAACCCTCAGCC
AAACAAGATGGACAGGCTGACGGAGATCCTGAACAGCATGAGGAATAACAGCAGTGACGTGGACGCCAAG
CTCACCTCCTTTCATGGAGGAAGCGCAGAACTCCACCAACTCCGAGGAGATGCTGGGAGAGATCGTGCGGA
CCATCTACCAGAAGGCCGTGTCAGACCGCAGCTTCGCCTTACGGCCGCCAAGCTCTGCGACAAGATGGC
CCTTTCATGGTCGAGGGGACCAAGTCCGGAGCCTGCTCCTCAACATGCTCCAGAAGGACTTCACTGTG
CGTGAGGAGCTACAGCAACAGGATGTGGAGCGCTGGCTGGGCTTACACCTTCTGTGTGAGGTGTTTG
GCACCATGCGGAGCAGCACTGGCGAGCCCTCCGAGTGTCTTGCAGTCTCAGGATGTGAAGGAAGA
CGCTGTCTCTGCTGCTCCATGGAGCTGCAAGCACAGGCCGGCTACTAGAGGAGCAGTCCCGGAGATG
ATGACGGAGCTCCTAGCCAGCGCGAGATAAGATGCTGTGCCCTCAGAGTCCATGCTAACACGCTCCC
TGCTCCTGGAGGTCATCGAGCTCCACGCCAATAGCTGGAACCCTCTGACGCCCCCATCACGCAGTACTA
CAACAGGACCATCCAGAACTGACAGCCTGA

ACGGTACGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001108891



[View online »](#)

Insert Size:	1011 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001108891.1</u> , <u>NP_001102361.1</u>
RefSeq Size:	1430 bp
RefSeq ORF:	1011 bp
Locus ID:	364900
Cytogenetics:	18q12.2