

Product datasheet for **RN208647**

Olr1585 (NM_001001109) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olr1585 (NM_001001109) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Olr1585
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN208647 representing NM_001001109 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGATCATGTCAATTGCACCTGGAACCTGGATCTTCATCCTTGCTGGTTTCAATACCTCTGGCACTCTGC
 AACAGCTTGCAGTCTTTGGGACCTGTGCATCTATCTCCTCACACTTGCAGGGAACCTGTTTCATCATTGT
 CTTGGTTTCAGGCAGATTCGGGGCTGTCCACTCCCATGACTTCTTTATCAGTGTCTCTCCTTCCTGGAA
 ATCTGGTATGTCAGCACCACAGTGCCACCTTGCTGCATACCTGCTCCACGGGCATTACCCATCCCTT
 CAGCTGCATGCTTCATCCAGCTGTATGCTTCCACTCCTTGGGCATGACTGAGTGCTACCTGTTAGGTGT
 CATGGCTCTGGACCGCTACCTTGCCATCTGTGCGCCATTGCACTACCATTCCTCTCATGAGCAGACAGGTA
 CAGAGAGAGCTAGTTCTGATTACTTGGGTGGCTGGCTTTTCAGCTGCCCTGGTGCCTGCAGTCTCACTG
 CCACTTTACCCTATTGTTTGAAAGAGGTAGCCCACTATTTCTGCGACCTGGGACCAGTGATGCAGTTGGC
 ATGTGTGGACACTAGCTGGCATGCTCGAGTTTATATTGCAGTGATTGGCATGATCAATACATGCAATCTC
 GTTTTCATCTTGGGACTCTATGGAGGTATTGTGAAAGCTGTGCTAAGGCTGCCTTCAGCTGCCAGTCGTG
 CCAAAGCTTTCTCCACATGTTCTCCACATCACTGTAGTGACACTGTTCTTTGGCTCTGCCTTCATTGT
 CTATGTAGGGCCACCAGAGATTCGAGCTGAGGGCAGAGGCAAGCTTATTGCTTTGGTATATACCTTGTTT
 ACCCTTTCTTCAACCCTATCATTATACTCTTCGCAATAAGGAAGTGAAGCGGCTTTCAAGAGGGTCA
 CACAGAGGATTAAGGGTGTGCTGGGT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001001109
Insert Size:	939 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001001109.1, NP_001001109.1</u>
RefSeq Size:	939 bp
RefSeq ORF:	939 bp
Locus ID:	406005
Cytogenetics:	13q24
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]