

Product datasheet for **RN208728**

Hapln4 (NM_001108398) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Hapln4 (NM_001108398) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Hapln4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN208728 representing NM_001108398 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGTGTGCTCCAGGAGCCCTCGGTACCGCGCACTCTGGGCCGTGGCCTGGGGACTCCTGCTGCTGG
TCCCGTGTCTCGCAGGGGCGCAGCGTGGCAGGAAGAAAGTCGTGCACGTGCTGGAGGGCGAGTCCGGTTC
GGTGGTGGTGCAGACGGCGCCGGGGCAGGTGGTGAGCCACCGCGGTGGCACCATTGTCTTACCCTGCCGT
TACCACACGAGGCGGGCGCCACGGTCACGACGGCGTTTCGCCTCAAGTGGACGAAAGTAGTAGACCCTT
TGGCTTTCGCGGATGTCTTCGTGGCTCTGGGTCCCCAGCATCGAGCCTTCGGTCCCTACCGAGGGCGGGC
CGAGCTACAGAACGATGGACCCGGGGACGCTCCCTGGTCTCCGAAATGTCACCCTGCAGGACTACGGG
CGCTACGAGTGCGAGGTACCAATGAGTTGGAGGATGACGCCGGCATGGTAAAGCTGGATCTAGAAGGGG
TGGTCTTCCCCTATCACCCGCGCGGAGGCCGCTACAAGATGACTTTCGCGGAAGCCAGCGCGCTTGTGC
GGAGCAGGACGGCATTCTGGCATCTGCGGAGCAGCTGCACACGGCCTGGCGCATGGCCTGGACTGGTGC
AATGCAGGGTGGCTACGCGATGGCTCCGTGCAGTACCCGGTGAGCCAGGCCCGAGAGGCCTGTGGCGGAT
CTGGGAGCCCTGGGGCTGGTGGCGGCACCAACGGGGCGTGCGCAACTATGGATATCGTCACAACGCCGA
GGAACGCTACGATGCCTTCTGCTTACGTCCAACCTCCCCGGCCGCGTGTTCCTGAAGCCGCTGCGA
CCTGTGGCCTTGGCCGGGGCCGTGCGTGTGCGGCACGGGTGCCACAGTGGCCAAGGTCGGGCAGC
TCTTCGCTGCGTGGAAGCTGCAGCTCTTGACCGCTGCACTGCGGGCTGGCTGGCGGATGGCAGCGCGCG
CTACCCCATCGTGAATCCGCGAACTCGCTGCGGAGGTCTCGGCCCGGCGTGCGCAGCCTCGGTTTCCCT
GACGCCTCCCGTCGCCTCTTCGGTGTCTACTGCTACCGCGCTCCGGGTGCCCCAGACCCAGCTCCAGGCG
GCTGGGGCTGGGGCTGGGCCGGCGGAGGAGTTGGGCCGAGGGTCGAGGGATCCGGCTGCCTGGACCCC
ACTGCGTGTCTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001108398
Insert Size:	1203 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001108398.1, NP_001101868.1</u>
RefSeq Size:	3687 bp
RefSeq ORF:	1203 bp
Locus ID:	361129
Cytogenetics:	16p14