

Product datasheet for **RN208735**

Calr3 (NM_001012212) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Calr3 (NM_001012212) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Calr3
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN208735 representing NM_001012212 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTCCGGGCTCGCGTGTGCTCTGGCCATCTGCGTGTGCGCGTGGCACTGGCCACTGTATACTTCC
AGGAGGAGTTTCTAGACGGAGAGCGCTGGAGAACCGCTGGGTGCAGTCCACCAATGACTCCCAGTTCGG
GCACCTCCGAGTCTCATCGGGAAAGTTCTATGGGCATAAAGAGAAAGACAAAGGCCTGCAGACCACCCAG
AACAGCCGATTCTACGCCATCTCTGCAAGCTTCAAGCCATTTCAGCAACAAGGGAAGACCCCTGGTAAATCC
AGTACACGGTGAAGCATGAGCAGAAGATCGACTGTGGCGGTGGCTACATCAAGGTCTTCCCGTCTGACCT
GGATCAGAAGAATCTGAATGAAAGTCCAATACTACATCATGTTTGGACCTGATATCTGTGGATTTGAT
ATCAAGAAAGTTCATGTTATCTTATATTCAAGAATCAATATCATGAAAAAAGAAGTCGATCAGATGCA
AGGTGGATGGCTTACGCATCTTACACCCTGATCTTAAGGCCAGATCTTCTTATGAGGTAAAAGTCGA
TGGCCAGTCCATCGAGTCTGGCAGCATTGAGTATGACTGGAACCTGACATCGCTGAAGAAAATGGAAAA
ACATCACTGGAGTCTAGGGACTGGGATGAGGTGGAAGGTAACAAAGCCCAGGACTGGGAAAAGCACTTTC
TGGATGCCGGTCCCAGCAAGCCAAGCGACTGGAACAGTGAGCTGGATGGGGACTGGCTGCAGAAGCCACC
ATACGAGGATGGCCTGAAAGCTGAGGGCATCGACAAGGACGTGTGGCTCCACCAGAAGATGAGGCCCGCT
GGTACCTGACGCAGTATGACCTCTCAGAGTTTGAGAATCGGGCCATTGGTCTGGAGCTGTGGCAGG
TGAGATCTGGAACCATCTTTGATAAATTCTGATCACAGATGATGAAGAATACGCAGAGAAGTTGGCAA
GGCCACCTGGGGGAAACTAAGGTTCCAGAAAAGGAGATGGATGCCATACAGGCCAAGGAAGAAGTAAAA
AAAGCCCAGAGGAAGATGAGGATCTGCTGATGGGGAAGTTTCACAGGCACAACCATTTTCAGCAGATTCC
ATAGGCAGGGCGAACTG**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001012212



Insert Size:	1140 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001012212.1</u> , <u>NP_001012212.1</u>
RefSeq Size:	1424 bp
RefSeq ORF:	1140 bp
Locus ID:	364529
Cytogenetics:	16p14