

## Product datasheet for **RN208747**

### **Raly (NM\_001011958) Rat Untagged Clone**

#### Product data:

|                           |                                                                                    |
|---------------------------|------------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                                |
| Product Name:             | Raly (NM_001011958) Rat Untagged Clone                                             |
| Tag:                      | Tag Free                                                                           |
| Symbol:                   | Raly                                                                               |
| Mammalian Cell Selection: | Neomycin                                                                           |
| Vector:                   | pCMV6-Entry (PS100001)                                                             |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                               |
| Fully Sequenced ORF:      | >RN208747 representing NM_001011958<br>Red=Cloning site Blue=ORF Orange=Stop codon |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTCCTGAAGATTCAGACCAGCAATGTAACCAACAAGAATGACCCTAAGTCCATCAACTCTCGGGTCT  
TCATTGGAACCTAAACACAGCTGTAGTGAAGAAGTCAGATGTGGAGACCATCTTTTCAAATATGGCCG  
TGTGGCTGGCTGTTCTGTGCACAAAGTTATGCTTTTGTCCAGTATGCCAATGAGCGCCATGCCCGAGCA  
GCTGTGCTGGGAGAGAATGGGCGGGTCTGGCTGGACAGACCCTGGACATCAACATGGCTGGAGAGCCCA  
AGCCTAATAGACCCAAGGGGCTAAAGAGAGCAGCAACTGCCATATACAGGCTGTTTGATTATCGAGGCCG  
ACTTTCTCCAGTGCCTGTGCCCGGGCAGTTCCTGTGAAGCGACCCCGTGTACAGTCCCTTTGGTTCGC  
CGTGTCAAACTACAATACCTGTCAAGCTCTTTGCCCGTTCCACAGCTATCTCTACTGGCTCAGCCAAGA  
TCAAGTTAAAGAGCAGTGAGCTACAGACCATCAAAACGGAGCTGACACAGATCAAGTCCAACATTGATGC  
CCTTTTGGGTCGCTTGGAGCAGATTGCTGAGGAACAGAAGGCCAATCCAGATGGCAAGAAGAAGGGTGAC  
AGCAGCAGTGGTGGAGGAGGGGGCAGCAGCAGTGGTGGAGGCGGCAGTGGCAATGTTGGTAGTGGCGGCA  
GCGGTGGCAGTGGGAGCAGCAGCAGCCGGCCACCAGCACCCCAAGAAGACACAGCTTCTGAGGCAGGCAC  
ACCCCAAGGAGAAGTCCAACTCGAGATGATGGTATGAGGAGGGACTGCTCACACATAGCGAGGAGGAG  
CTGGAGCACAGCCAGGACACAGATGCAGAAGATGGGCCTTGCA**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

|                    |              |
|--------------------|--------------|
| Restriction Sites: | SgfI-MluI    |
| ACCN:              | NM_001011958 |
| Insert Size:       | 888 bp       |



[View online »](#)

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).                                                                                                                                                                    |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                              |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <u>NM_001011958.1, NP_001011958.1</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>RefSeq Size:</b>           | 1509 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>RefSeq ORF:</b>            | 888 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Locus ID:</b>              | 296301                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Cytogenetics:</b>          | 3q41                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |