

Product datasheet for **RN208819**

Neil1 (NM_001025754) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Neil1 (NM_001025754) Rat Untagged Clone
Tag: Tag Free
Symbol: Neil1
Synonyms: MGC116267
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN208819 representing NM_001025754
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCCGGAGGGCCAGAGCTGCACCTGGCCAGCCACTTTGTGAATGAGATATGTAAGGGCTGGTATTCCG
GTGGGTGTGTGGAGAAGTCTCTGTGAGCCGAAACCCGGAGGTGCCCTTTGAGAGCAGTGCCTACACAT
CTCAGCTTTAGCCCGAGGCAAGGAGCTGCGCTTGACATTGAGCCCTGCCTGGGTCCCATCCCCCTCAG
AAACCACTGTCCCTTGTCTCCGCTTTGGTATGTCTGGATCCTTCCAGCTGGTACCCGCAGAAGCACTGC
CGCGCCACGCCATCTACGTTTCTATACGGCCCACTGCCCGCCGGCTTGCCTTTGCTTCGTAGACAT
CCGCCGCTTTGGCCACTGGGACCCTGGAGGTGAATGGCAGCCAGGCCGCGGACCCTGTGTCTTGTGGAG
TATGAACGGTTCAGAGAGAACGTGCTTCGGAACCTATCGGACAAAGCCTTTGACCGGCCATCTGCGAGG
CCTTGTGGACCAGAGTTCTTCAATGGAATCGGCAACTATCTGCGGGCAGAGATCTTGTATCGGCTGAA
GATCCCCCTTTTGAGAAGGCCGTACAGTTCTAGAGGCCCTGCAACAGTGCCGGCTGAGCCAGAGCTG
ACCCTGAGCCAGAAGGTCAAGGCCAACTGCAGAACCCAGACCTGCTGGAAGTGTGCTGCTGGTCCCA
AGGAAGTGGTTCAGCTGGGGGGCAAAGGCTATGGCCAGAGCTGGAGAGGAAGATTTTGTGCTCCCTTCG
AGCCTGGCTTCGGTGCTATGCTGTGCCAGGCATGAGCTCCCTGAGAGACCGGCATGGCCGTACCATTTGG
TTCCAGGGTGATCCTGGACCCTTGGCACCCAAAGGGGCGAGATCCCGAAAAAAGAAGTCACAGGAAACAC
AGCTGGAGGCTGAGGACAGAAAAGAGGACCCTCCAAGCAAGTCCATTTCCAGGACTCAGAGGGCCAGGAA
GCACTGTCTAAGAATGTACCTCAACAGTCTGAAGGGGCCAGCCTCCACAAAACCTGGAAACCCCTGCA
GTCCCTGAGAAAGGAAAGAGAAGGGGGAGACGGGCGAGCACAGGCCGCCACAGAGGACCAAGAGTAAGC
CTGATACTCAACCCAGGGAGGCTGGGGAGATTTAGCCTCATAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online >](#)

ACCN:	NM_001025754
Insert Size:	1164 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001025754.1, NP_001020925.1</u>
RefSeq Size:	1768 bp
RefSeq ORF:	1164 bp
Locus ID:	367090
Cytogenetics:	8q24