

Product datasheet for **RN208843**

Actrt3 (NM_001108942) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Actrt3 (NM_001108942) Rat Untagged Clone
Tag: Tag Free
Symbol: Actrt3
Synonyms: Arpm1; RGD1561457
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN208843 representing NM_001108942
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAGCAGCTACCAGCCACCAGTGGTATCGACAACGGCTCGGGAATGATCAAGGCGGGCCTGGCTGGGG
 CCCGGGAGCCCCAGTTCGTCTACCCGAATATTCTGGGCCGGACCAAGAACCACAGCCCTGCGGCTGACAG
 CAAGCAGGAGCTGCGTGTGGGTGATCAAGCGCAGGAGCGCAGGAGCTTCTGTCCATCAGCTATCCAGTG
 GAACGGGGTCTCATTTCCTCTGGGAGACATGGAGATCATGTGAAACATATCTATGACTACAACCTGA
 ACCTGAACCCAGCGATGGCCCGTCTTGTTACAGAACCGGCTCAACCCATTGGCAGACCGGCAGCA
 CATCTCCGAAGTGTTCCTGAAAATCTGGGAGTCCCTGCCTTCTATATGTCTGTCCAGGCTGTGTTGGCT
 CTGTTTGCTGCTGGTTTTACCACTGGCCTGTGCTGAACTCAGGTGCTGGGATTACTCAGTGTGTGCCCA
 TCTTTGAGGGTTACTGCCTATCACATGGTGTAAAAACAGCTAAACGTGGCAGGAATTGACCTCACCAACTA
 CCTCATGATGTTGCTGAAGGATGATGGTATCATGCTGCTTAGAACTGGGGACAGAAAACTGTCACAGAC
 ATTAAGGAGAATGCTTGCTATGTGGCAATGAACTATGATGAGGAAATGGTCAAGGAATCTAACGTAGAGA
 AAATGTACACCCTCCCTGACGGGAAGACTGTAAACTCCGTAAGCAGGTATTCGCTGCCCGGAGGCCCT
 CTTCTCCTCGTATCTTGTAATGTGGATGCTCCTGGCATCGACAGGATATGCTTCAGCAGTATAATGAAG
 TGTGATGCAGATCTCAGGAACCTCTTCTCCAACATCATCCTTCCGGAGGATCGACCTCTTTTCTCTG
 GTTTAGATAAGCGACTCATTAAAGGATGTGGCAAAGTTGGCAACCTGCCAATACCACCGTGCAGGTTGTAGC
 TCCTCCTGAAAGGAAAATTTTCGGTGTGGATGGGGGATCTATTCTGTCATCCTTGTCTGCCTTCCAAGAC
 ATGTGGATCACGGCTGCAGAATTCGAAGAAGTGGGACCAACATTGTACACCAGAGGTGCTT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online >](#)

ACCN:	NM_001108942
Insert Size:	1116 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001108942.1</u> , <u>NP_001102412.1</u>
RefSeq Size:	1116 bp
RefSeq ORF:	1116 bp
Locus ID:	365763
Cytogenetics:	2q24