

Product datasheet for **RN208854**

Ankrd37 (NM_001108400) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Ankrd37 (NM_001108400) Rat Untagged Clone

Tag: Tag Free

Symbol: Ankrd37

Synonyms: RGD1565715

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN208854 representing NM_001108400
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCTGTTGCTTAGCTGCAACTTGGAGGAGGATGGCCTTAAGAATTTGCTGGAGACAGGAGCATCAGTCA
 ACGCACCCCGGATCCTCAAGAGCAGTCGCCTGTTACCTAGCTGCAGGTGGCGCCTTGCTTGCTTTCT
 TCTCTGGCAGCTGCAAACCGGTGCTGACCTCAACCAACAGGATGTGTTGGGAGAACTCCACTCCACAAG
 GCAGCCAAAGTTGGAAGCCTAGAGTGCCTTAGCCTCCTGGTGGCCAGTGATGTCCGGATTGACTTATGTA
 ATAAAAATGGGCAAACAGCTGAAGATCTTGCTTGGTCATGTGGCTTCCAGAATGTGCCAGGTTTCTAAC
 AATGATTAATGTACGCAGACAGCAAAATCCTCACGTGAGCACCTGAGGGTGACCCCTATGCTCCTGTG
 CTGAGACAGAAGCGGAGTTTCAGAACTATGGAGAGCAGAGTGTTGAAGAGAAAGTGT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001108400

Insert Size: 480 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001108400.2, NP_001101870.1</u>
RefSeq Size:	793 bp
RefSeq ORF:	480 bp
Locus ID:	361149
Cytogenetics:	16q11