

Product datasheet for **RN208926**

Lrp11 (NM_001106217) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Lrp11 (NM_001106217) Rat Untagged Clone
Tag: Tag Free
Symbol: Lrp11
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN208926 representing NM_001106217
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGGTGCCTCAGCCAGGAACCCTGAAGCTGTCCACCTGAAGGAAGGAGCCTACGTCTTCCAGCTGA
CTGTGACTGACTCCGTGGGTCAGAGAAGCTCCGACAATGTATCTGTGACTGTGCTTCCAGGCTCCGCTC
CACGGGAGGATGCTCCAGTGCTTGTCCGATACCATTTCTTCTGTGACAATGGCTGCTGCATTGACATC
GCCTTGGCGTGTGATGGCGTGAGGCAGTGTCCAGACGGTTCTGACGAAGACTTCTGTCAAACCTGGCTA
CTTCTGCCAGCCAGGAGCCACAGGATTAATGAAGGTGAAGGAGACCCCTGGTTGGAGAAATCCCAGAA
AGCCACTACCCACAACCAACCAGCTACTGTATCCACCCAGAGAAGAGGATCCAGCCACCCAGAAGGCCG
CCAGAGAGTCAGATCAACCCTGTCCAGCCAGACAGTAATTCCTCAGGAAAAGAAATCAAGAAGAAAAATACG
ATCTCAAGTCTAAGAGTGACCAAGCACGAGGGGAGCACCCAGCCAGAGCAGGTGCAGTGTGCCCT
GGCACTGGGTCTGGCCATCACTGCTTTGTTACTTCTCATGGTGACATGCCGACTCCGACTGGTGAAGCAG
AAACTAAGAAAGCCGTCCTTACATCTGAGGAATCCGACTACCTCATCAATGGGATGTATCT**GTA**A

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja1338_e07.zip
Restriction Sites: SgfI-MluI
ACCN: NM_001106217
Insert Size: 699 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001106217.1</u> , <u>NP_001099687.1</u>
RefSeq Size:	3255 bp
RefSeq ORF:	699 bp
Locus ID:	292462
Cytogenetics:	1p13