

Product datasheet for **RN208939**

St6galnac6 (NM_001015036) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: St6galnac6 (NM_001015036) Rat Untagged Clone
Tag: Tag Free
Symbol: St6galnac6
Synonyms: Siat7F
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN208939 representing NM_001015036
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAACGATGCTCCACCACTGGCTACTCGGCCGACGTCGGCAACAAAACCACCTCCGCGTAGTGGCTC
ATCCAGTGTGTCCGGGTGCTGCGGAAGCCCCAGGAATTTGTCAACAGGACCCCTGAAACTGTGTTTCAT
CTTCTGGGGCCCCCAACAAGATGCAGAAGCCCCAGGGCAGCCTCCTGCGGTCATCCAGCGGGCGGGC
CTTGCGTTCCTAACATGGAGGCCATGCGGCTCTCAAGCCCGCATGCAGCAGTTTGACGACCTCTCC
GGGGTGAGACAGGCAAGGACAGGAAAAGTCCCACTCGTGTTGAGCACAGGCTGGTTTACCATGGTGAT
TGCGGTGGAGTTGTGACCATGTGCACGTGTATGGCATGGTCCCTCCGACTACTGCAGCCAGCGGCC
CGCCTGCAGCGGATGCCGTACCACTACTATGAACCAAGGGGCCAGATGAATGTGTACCTACATCCAGA
ACGAGCACAGCCGTAAGGGCAATCACCACCGCTTCATCACTGAGAAGAGGGTCTTCTCATCTGGGCTCA
ACTCTACGGTATCACCTTCTCCACCCCTCCTGGAC**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001015036
Insert Size: 600 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001015036.1</u> , <u>NP_001015036.1</u>
RefSeq Size:	2243 bp
RefSeq ORF:	600 bp
Locus ID:	407765
Cytogenetics:	3p11