

## Product datasheet for **RN208945**

### Serpinb6b (NM\_001012214) Rat Untagged Clone

#### Product data:

|                           |                                                                                    |
|---------------------------|------------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                                |
| Product Name:             | Serpinb6b (NM_001012214) Rat Untagged Clone                                        |
| Tag:                      | Tag Free                                                                           |
| Symbol:                   | Serpinb6b                                                                          |
| Mammalian Cell Selection: | Neomycin                                                                           |
| Vector:                   | pCMV6-Entry (PS100001)                                                             |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                               |
| Fully Sequenced ORF:      | >RN208945 representing NM_001012214<br>Red=Cloning site Blue=ORF Orange=Stop codon |

TTTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGGATCCTCTGCTGGAAGCAAATGGCACATTTGCCTTCAACCTTTTGAAAACACTGGGTGAAGACAGCT  
CAAAAAATGTACTTTTCTCACCCTGAGCATCTCCTCAGGCCTGGCCATGGTCTTCATGGGGGCAAAAGG  
AACAACTGCACATCAGATGATTCAGGCACTTTCTTTGGATAAATGCAGCGGCAGAGGAAGTAGAGATGTT  
CACCAGGGCTTCCAGTCCCTTCTCGCCAAAGTGAACAAGACTGGCACACAGTACTTGCTCAAAACAGCCA  
ACAGGCTCTTCGGGGAAAAGACTTTTGATATCTTAGCCTCTTTAAAGATGCCTGCCGCAAGTTCTACGA  
AGCAGAAATGGAAGAGCTGGACTTCAAGGGTGCTCCAGAGCAGTCGCGACAGCACATCAACACCTGGGTA  
GCCAAAAAGACAGAAGAAAAAATTACCGAGTTGCTATCTTCAGGTTCAAGTGAATGCAAACTCCACTGG  
TCCTTGTGAATGCCATCTACTTCAAAGGGAAGTGAAGAAGCAGTTTAAACAAAGAGGACACTCAAGAAAT  
GCCTTTCAAAGTCACCAAGAATGAGGAGAAACCTGTGAAATGATGTTTAAAGAACTACCTTTAAGATG  
ACCTATGTGGAAGAGATATCCACTACTATTCTGTTGCTTCCCTATGTTGGCAATGAACGAACATGATCA  
TCATGCTCCAGATGAACACATTGAAGTGAAGTGGGAAAAGGAAATAACTTACAAGAAATTCATAGA  
GTGGACGAGCCTGGACAAGATGGAAGAAAGAGAGGTGGAGTTTTCTCCCAAAGTTTAACTGGAGGAG  
AACCATGACATGAAGGATGTCCTCCACAGGCTGGGCATGACTGATGCCTTTGAGCAGGGCATGGCAGACT  
TTTCAGGAATAGCATCTAAAGAAGGCTTGTCTGTCCAAGGTCAATCACAAGTCCTTTGTGGAGGTCAA  
TGAAGAGGGCACAGAGGCTGCAGCTGCCACTGCAGCCAATGTAACTTTTAGGTGTATGGTCCCTTACTTT  
TGTGCCAACCAACCTTTTCTTTCTTATCCAGCACAGCAGGACCAATGGGATTGTATTCTGTGGCCGGT  
TCTCCTCCCCA**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>ACCN:</b>                  | NM_001012214                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Insert Size:</b>           | 1134 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).                                                                                                                                                                          |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>RefSeq:</b>                | <u>NM_001012214.1, NP_001012214.1</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>RefSeq Size:</b>           | 1436 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 1134 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>              | 364705                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Cytogenetics:</b>          | 17p12                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |