

Product datasheet for **RN209025**

Ier5l (NM_001025041) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ier5l (NM_001025041) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Ier5l
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN209025 representing NM_001025041 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAGTGCGCCCTGGACGCCAGAGCCTGATCAGCATCTCCCTGCGCAAGATCCACAGCTCCCGGACCC
AGCGCGGCGGCATCAAGCTGCACAAGAACCTCCTGGTGTCTACGTGCTCCGCAACGCGCGCCAGCTCTA
CCTGAGCGAGCGCTACGCCGAGCTCTACCGCGCCAGCAGCAGCAGCAGCAGCAACAGCAGCCGCC
CACCACCAGCACCAGCACCTCGTTACGCGGCGCCCGGAATGCCGGCCAGCGCGGCCGACTTCGGCCCGC
TCCAACCTGGTGGCGGCGGGACGCGGAGGCGCGAGCCGGTCGCCCGGCACCAGCTGCACCAGCTCCA
CCAGCTCCACCAGCTGCACCTCCAGCAGCAGCTGCACCAACACCAGCACCCGGCGCCAGGGGCTGCACG
GCGGCGGCGCCGGTGGCGGTGGCCGGGCGCCCGCGGGCTGCGCGGGGGCGCTCTCGGAGCTGCCCGGGT
GCGCTGCGCTCCAGCGCCGCACGGCGCGCCCCACCGCGGGCAGCACTTGGAGCCGCTGCAGCCCGGTCC
TGCGCCGCTGCCGCCGCCCGCGCCCGCGCTCTGCCCGGGGACCTCGCGTCCCGGCCGCTGCTCT
GCGCCCTCGGGCGCTCCGGGGGCTGCTCTCGACCGTCGCTGCCTCCTCTCGGCTGCCTCTCCAGCCC
CTTCTCATCCCCTGGCTTCTACCGGGGCGGTACCGGCCCCCTCGGACTTCGGCGTACACTGCAGCAG
CCAGACCACTGTGCTGGACCTGGACACTCACGTGGTGACCACGGTAGAAAACGGCTACTGCACCAGGAC
TGCTGCGCCTCCGCCACTGCCCTGCTGTGGCCAGGGCGCTCCGGGACCGGGTCTGGCGTCCGCCGCTG
GTTGCAAGCGCAAGTATTACCCTGGCCAGGAGGAAGAGGACGACGAGGAGGAGGATGCGGGCGACCTGGG
AGCCGAGCCCCCGGGGGCAGCCGTTGCCCCCTGTAAGCGCGCCCGCTTCGAGGACTTCTGCCCGGAC
TCGTCCCCGGACGCGTCCAACATCTCAAACCTGATCTCCATCTTTGGCTCGGGCTTCTCGGGGCTGGTGA
GCCGACAGCCGGACTCCTCGGAGCAGCCGCCCGCTCAACGGGCAGCTGTGCGCCAAGCAGGCGCTCGC
CAGCCTCGGCGCTGGAATCGAGCCATTGTCGCTTC**TAG**

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-RsrII



[View online »](#)

ACCN:	NM_001025041
Insert Size:	1230 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001025041.1, NP_001020212.1</u>
RefSeq Size:	1612 bp
RefSeq ORF:	1230 bp
Locus ID:	499772
UniProt ID:	<u>Q5PQP0</u>
Cytogenetics:	3p12