

## Product datasheet for RN209055

### Ebi3 (NM\_001109421) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ebi3 (NM\_001109421) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Ebi3  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN209055 representing NM\_001109421  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGCATCGCC

ATGTCCAAGCTGCTCCTCTGTCACTTGTCTCTGGCCAGCTGCTCCCCTGGTAACACTGGGACAGCGC  
 TAGTGGATCTGAGCCAGCCAGAGTGCATGCCATGCTTCTCGGTACCCGTGGCTGTGGACTGCTCCTG  
 GATGCCACTCCAGGCTCCCAACTCCACCAGAGCCATGTCCTTCATTGCCACTTACAGGCTCGGTGTGGCC  
 GCCCAGCAGCAGAGCCAGCCCTGCCTGCAACGGACCCCCCAGGCATCCCGTTGCACCATCCCCAGGTGC  
 ACCTGTTTTCCACAGTGCCTACTTGTAAATGTCACTGTAGTGACCCCTGGCGCGCCAACAGCAGCCT  
 TCTAGCCTTTGTGGCGGAACGAATCATCAAGCCGACCCTCCGGAAGGCGTGCGCCTGCGCAAAGCGGGA  
 CAGCGCTGCAGGTTCTCTGGCATCCCCCTGCTTCTGGCCCTTCCCGGACATCTTCTCTCAAGTACC  
 GACTTCGCTACCGGCGCCGAGGAGCCTCTCACTCCGCCAGGTGGGACCTATTGAAGCCACGACTTTCAC  
 CCTCAGGAACACACAACCCCATGCCAAGTATTGCATCCAGGTGTCAGCCAGGACCTCACCGATTATGGG  
 AAACCAAGTGACTGGAGCCTCCCTGGGCAAGTAGAGAGGGCGCCCCAGAAGCCCTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001109421  
**Insert Size:** 687 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001109421.1</a></u> , <u><a href="#">NP_001102891.1</a></u>
<b>RefSeq Size:</b>	1198 bp
<b>RefSeq ORF:</b>	687 bp
<b>Locus ID:</b>	680609
<b>Cytogenetics:</b>	9q12