

## Product datasheet for **RN209098**

### Cited1 (NM\_172055) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Cited1 (NM\_172055) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Cited1  
**Synonyms:** Msg1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN209098 representing NM\_172055  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGCCAATATGTCGAGGCTGCACTGGATGTCAAGGGTGGTACCACCTCTGAAAGGAAGATGCCAAC  
AGGAGATGAACTCTTGGCTACTCTAACCTCGGGTAAAAGATCGCAAGGCAGTGACCGTCTGCACTA  
CCCCGGGGTACCGCAAATGGAGCAAAGCCAATGGAGTGCCACTAGCCCTCTGGATCATCATCTCCA  
ACAGGCTCTCTACTGCCACCCCTCTACCAAACCCCATCCTTCAACCTGCACCCACCCCTCACCTGC  
TGCCAGTATGCAGCTTCAAGCTTAATAGCCAGTACCAAGGGGCTGCAGCGACTGCTGCTGCTGCCCT  
GGCTGCACCAAGCCAACCAGGAGAGGATGAGCCCTGCTAACTGGGGCACTGGGATCCAGGCAGGAGCA  
GGGGATCACCCGGATCAGTCTCTCTGCTGGTGGCCAGAGCCCTGCTATCATTGATTCTGACCCAGTGG  
ATGAGGAGGTGCTGATGTCTCTGGTGTAGAATTGGGGCTAGACCGAGCCAATGAGCTTCCCGAGCTGTG  
GCTAGGGCAGAATGAGTTTGATTTCACTGCAGATTTTCCCTCTGGCTGC**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_172055  
**Insert Size:** 612 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_172055.2</a></u> , <u><a href="#">NP_742052.2</a></u>
<b>RefSeq Size:</b>	873 bp
<b>RefSeq ORF:</b>	612 bp
<b>Locus ID:</b>	64466
<b>Cytogenetics:</b>	Xq22
<b>Gene Summary:</b>	mouse homolog acts as a transcription factor; may interact with p300/CBP, TATA-binding protein, and Lhx2 [RGD, Feb 2006]