

Product datasheet for **RN209104**

Fra10ac1 (NM_001014246) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fra10ac1 (NM_001014246) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Fra10ac1
Synonyms:	RGD1309482
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN209104 representing NM_001014246 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCATGGGCATGGGGCTATGACTCTGATTTTAGTGATGATGAGCAAGGTGGAGGATCCAGCAAAAAGA
GGAAGAAGACAGTTGAAGATGAGTTGCTGCTTACAAAACCATTTCAAAAAGAAAAATATGGAAAGGTGGC
CCATAAACAGTTGCTGCAGACTTGCTGGACAGGAAGAAGCAAGAAATAGAAGGTTTCATCTCATAGCT
ATGGATGCTTATCAAAGGCACACAAAGTTGTAATGACTATATTTGTACTATGGTGGCAAAAGAGAGG
ACTTTAAGCGGTTGGGGAGAATGATAAGACAGATCTGGATGTTATACGAGAAAATCACAGATTCCTATG
GAATGAGGAGGATGAAGCAGACATGACTTGGGAGAAGAGACTTGCAAAGAAATACTATGATAAATTATTC
AAGGAATATTGCATAGCAGATCTCAGTAGATACAAAGAAAATAAGTTTGGATTTAGTGCGAATAGAAA
AAGAAGTAATTTCAAGAAAAGGCCAGTTTTTTTGTGAAAATAAGTGTGATGAAAAAGAAGGTTTGAA
GAGCTGGGAAGTCAACTTTGGTTATACTGAGCATGGTGAAGAGAAATGCGCTTGTTAAGTTAAGATTA
TGCCAAGAATGTTCCTTAAATTTTATCACAGGAGAAAAGAAATCAAGTCCACAAAAAAAAGGT
CTAAAACCAAGACTGAGTCTGACGAGTCACCTCATAAGAACTCCAGATCCTCATCTTCAGAGGAGGCCTC
CCAGGGGAAGGATGAAGGACATTCATCCTCGAAAAGGTCAGAAGATTCTAGAAACAGAAACGCTGGTGAG
GAAGACAGTGCTTCAGACTCTGAGCTCTGGAAGGTCCTGCCCAGAGACAGATGAAAAGTCACAGGAAG
AAGAATTTGATGATTATTTTCAGGATTTGTTCTG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001014246



[View online »](#)

Insert Size:	948 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001014246.1</u> , <u>NP_001014268.1</u>
RefSeq Size:	1388 bp
RefSeq ORF:	948 bp
Locus ID:	365458
UniProt ID:	<u>Q5FVF1</u>
Cytogenetics:	1q53