

## Product datasheet for RN209140

### Kcnip4 (NM\_181365) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Kcnip4 (NM\_181365) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Kcnip4  
**Synonyms:** Kchip4  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN209140 representing NM\_181365  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAACGTGAGAAGGGTGGAAAGCATCTCGGCTCAGCTGGAGGAGCGAGCTCCACAGGCGGTTTCCTCT  
ACGCTCAGAACAACACCAAGCGCAGCATTAAAGAGCGGCTCATGAAGCTTTGCCCTGCTCAGCTGCCAA  
AACGTCGTCTCCTGCTATTCAAACAGTGTGGAAGATGAACTGGAGATGGCCACTGTCAGGCACCGGCCT  
GAAGCCCTGGAGCTGCTGGAGGCCAGAGCAAATTCACCAAGAAAAGAGCTTCAGATCCTTTACAGAGGAT  
TTAAGAATGAATGCCCCAGTGGTGTGTTAATGAAGAAACCTTCAAGGAGATTTACTCGCAGTTCTTCCC  
ACAGGGAGACTCCACCACATATGCACATTTTCTTCAATGCATTCGACACGGACCACAATGGAGCTGTG  
AGCTTTGAGGATTTTCATCAAAGGTCTTTCCATTTTGTCTTCGAGGGACAGTACAAGAAAACTGAACTGGG  
CATTTAATTTGTATGACATAAACAAAGATGGCTACATCACTAAAGAGGAAATGCTGGACATAATGAAAGC  
AATTTACGACATGATGGGAAAATGCACATACCCTGTCTCAAGGAAGATGCTCCCGACAGCACGTGGAG  
ACATTTTCCAGAAGATGGACAAGAATAAAGATGGTGTGTTACCATAGACGAGTTCATTGAAAGTTGCC  
AAAAAGATGAAAACATAATGCGCTCCATGCAGCTCTTTGAAAATGTGATCTAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_181365  
**Insert Size:** 753 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_181365.2</a> , <a href="#">NP_852030.1</a>
<b>RefSeq Size:</b>	1289 bp
<b>RefSeq ORF:</b>	753 bp
<b>Locus ID:</b>	259243
<b>UniProt ID:</b>	<a href="#">Q99MG9</a>
<b>Cytogenetics:</b>	14q11
<b>Gene Summary:</b>	a potassium channel interacting protein which modulates the opening and closing, and disrupts rapid inactivation of the Kv4 potassium channel [RGD, Feb 2006]