

Product datasheet for **RN209233**

Pttg1ip (NM_001013238) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Pttg1ip (NM_001013238) Rat Untagged Clone
Tag: Tag Free
Symbol: Pttg1ip
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN209233 representing NM_001013238
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCGTCTGCTGTCTTGGGCTGACGCTGCGCTGGGTGATGTTCCCTCAGTGCCGTGCTGCTTGTCTTC
 TGCCCGGAGCCTCTGCGCAGGAGCCTCCAGGAGTGGGTGCTCTGAGTACACAAACAGATCCTGTGAAGA
 GTGCCTCAGGAATGTCTCGTGTCTCTGGTGCAACGAGAACAAGCGTGCCTGGACTACCCAGTGAGGAAA
 ATCTTGCCCTGCTCTCTTTGTAATTGAGCTCCGCTCGCTGGGGCGTATGCTGGGTGAACTTTGAGG
 CCCTGATCATCACCATGTCGGTGTCTTGGGGCTCTGTGCTCTGGGCATCACCCTGTGCTGCTGCTG
 CTGCCCGGAAAAGAGAAGCCGGAAGCCAGACAAGAGCGATGAGCGGGCCATGAGAGAGCAGGAGGAGG
 AGAGTGCACAGGAGGAGAGGAGAGCAGAGATGAAGTCAAGACACGATGAAATCAGGAAGAAATATGGTC
 TGTTAAAGAACAAGAACCATATGAGAAGTT**CTGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001013238
Insert Size: 525 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001013238.1, NP_001013256.1</u>
RefSeq Size:	2259 bp
RefSeq ORF:	525 bp
Locus ID:	365548
UniProt ID:	<u>Q6P767</u>
Cytogenetics:	20p12
Gene Summary:	May facilitate PTTG1 nuclear translocation.[UniProtKB/Swiss-Prot Function]