

## Product datasheet for **RN209316**

### **Tm4sf5 (NM\_001107015) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids

**Tag:** Tag Free

**Symbol:** Tm4sf5

**Mammalian Cell Selection:** Neomycin

**Vector:** pCMV6-Entry (PSI00001)

**E. coli Selection:** Kanamycin (25 ug/mL)

**Fully Sequenced ORF:** >RN209316 representing NM\_001107015  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGC**

ATGTGTACTGGAAAGTGTGCTCGATGCCTGGGGCTCTCCCTCATCCCTCTCTCCCTGATCTGCATCGTGG  
CCAATGCCCTCCTGCTGGTACCTGATGGGAAGACCACCTGGACGGATGGCAACCTCAGCTTGCAAGTTTG  
GCTCATGGGCGGCTTCATTGGAGGGGGCTTGATGGTGTGTGCCCTGGAATTGCTGCAGTCCGGGCGGGG  
GGAAAGGGCTGCTGCGGTGCAGGCTGCTGTGGGAACCGCTGCAGGATGCTGCGCTCTGTCTCTCTCCG  
CCTTTGGGATACTTGGTGTATCTACTGCCTGTGAGTGGCAGGAGCTGGGCTCCGAGTGGGACCAATG  
CTTAATAAAATAACAAATGGGACTACCACTTCCAAGAAACCAAGGCGTTACTTGACAATGACTCTT  
TGGAATTTATGTGAGGCGCCACCTCACGTGGTACCCTGGAATGTGACGCTCTTCTCAATTCTGGTGGTGC  
CCTCAAGTCTGGAAACAGTGTGTGCGGAATACAGCTGGTGAATGCGACCTTGGGCGTGTGTGTGGCGA  
TTGCCGAAAAAGGAGGGCTCAGCTCACT**GA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_001107015

**Insert Size:** 591 bp



|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).                                                                                                                                                                          |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>Note:</b>                  | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>RefSeq:</b>                | <u><a href="#">NM_001107015.1</a></u> , <u><a href="#">NP_001100485.1</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>RefSeq Size:</b>           | 951 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>RefSeq ORF:</b>            | 591 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Locus ID:</b>              | 303256                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Cytogenetics:</b>          | 10q24                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |