

Product datasheet for **RN209329**

Cdk5r2 (NM_001109309) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Cdk5r2 (NM_001109309) Rat Untagged Clone

Tag: Tag Free

Symbol: Cdk5r2

Synonyms: RGD1562530

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN209329 representing NM_001109309
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGGCACGGTGTCTGCTTTCCCTGCTTCCTCGGCCAAGGGCCGGAGGCCGGGAGGGCTGCCGGAGG
AGAAAAAGAAAGGCGCGCCCGGGGGACGAGGCGCTGGGGGGCTACGGCGCGCCCGCGCCGCAAGGG
CGGCAAGGCGAGAGCCGGCTCAAACGGCCGTCCTACTCATCTCGGCGCTCACCTGGAAGCGCCTAGTG
GCAGCCTCGGCCAAGAAGAAGAAAGGCAGCAAAAAGGTGACGCCAAGCCAGCGTCCACCGGTCCCGACC
CTTTGGTCCAGCAACGCAACCGGGAGAACCTTCTCCGCAAGGGTGCAGACGCCCCGACGGCGGGGAGC
AGCCAAGCCCCTGGCCGTGCCGGTGCCACGGTGCCCGCTCCGCTGCCACCTGCAGCCCCCGTGGGG
GGCAGCGCGGCCGCCGCCACAGGCTCAGGCGGGGAAAGCCACCTCCGCCGCCACCTCCAGCCCCGC
AGGCGGCACCGCCAGCGCCCGAGGCTCGCCGCGACGGGTATCGTGCAGGCGTCCACCGCGAGCTGCT
CCGCTGCCTCGGCGACTTCGTGTGCCGACGCTGTACCGCTCAAGGAGCTGAGTCCCGGCGAGCTGGT
GGCTGGTTCGCGGCGTGGACCGCTCGTGTGCTGTCAGGGTGGCAAGACCAGGCCTTATTACCCCCG
CCAACCTGGTGTTCGTGTACCTGCTGTGCCGCGAGTCTGCGCGGGGACGAGCTGGCGTCCGCCGCCGA
ACTGCAGGCCGCTTCTCACCTGCCTTACCTCGCCTACTCTACATGGGCAACGAGATCTCTACCCA
CTCAAGCCCTTCTCGTGGAGCCGACAAGGAGCGCTTCTGGCAGCGCTGCCTGCGCCTCATCCAGCGGC
TCAGCCCGCAGATGCTGCGGCTCAACGCCGACCCCACTTCTTACGCAAGTCTTTCAAGACCTCAAGAA
CGAGGGCGAGGCCGCGGCGAGCACCGGGGTCCACCGAGCGGAAGCGGTCCACCACTCTCGTCGACC
GCCAGGGACAGCTGCGCGCCGGAGCCAAGCATTGGACTATGAATTGGACCGCTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001109309
Insert Size:	1107 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001109309.1, NP_001102779.1</u>
RefSeq Size:	1469 bp
RefSeq ORF:	1107 bp
Locus ID:	501164
Cytogenetics:	9q33