

## Product datasheet for **RN209365**

### Abi3 (NM\_001013118) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Abi3 (NM\_001013118) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Abi3  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN209365 representing NM\_001013118  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCGGAGCTGCAGCAGCTGCAGCAGTTGCAGGAGTTCGATATCCCCACGGGCCGGGAGGCTCTGCGGG  
GCAACCACAGCGCCCTGCTACGGGTGGCCAACTACTGTGAGGATAACTACTTGCAGGCCACAGACAAGCG  
GAAGCGCTGGAAGAGACGATGGCTTTCACCACCCAGGCCCTGGCCAGTGTAGCCTATCAAGTGGGTAAC  
CTGGCGGGGCACACGCTTCAATGCTGGATCTACAGGGTCTGCCCTGCGGCAGGTGGAAGCCAAGATGA  
GCACACTGGGCCAGATGGTGAACATGCACCTGGAGAAAGTAGCCAGAAGGGAGATTGGCACGTTGGCCAC  
TGTCTGCGGCTGCCCTAGCCAGAAGGTCATCCCTCCTGAGAGCCTGCCTCCCCTCACTCCCTACTGC  
AGAAAACCCCTCAACTTTGCCTGCTGGATGATGTTGGCCATGGAGTCAAGGACTTGAGCACACAGCTGT  
CACGGACCGGGACCCTGTCTCGAAGAGCATAAAGGCGCCGCTACACCTGCCTCTGCCACGCTGGGGAG  
ACCACCCCGGATCCCTGAGCCGGTGCAGCTCCCAGCGGTGCCAGACGGCAAGCTCTCCGCTGCCTCCTCT  
GTGTCTTCTTGGCCTCCGACGGCAGTGCAGAAGGTGCCAGTGGGATCCCCAGTCCAAGGGACAGGTAG  
CACCTGCAACCCCGCCTCCTCCACCTATAGCGCCTGTAACCTCCACCTCCTCCACCATTGCCTGCTGAGAT  
CTTCTTGCTGCCCTCCGATGGAGGAGTCCCAGCCCCCTCCGGAACAGAGTTGCCCTGCCTCCTCCT  
CCGGCTCTACAGGGGGATGAAGTGGGGTCTGCCTCCGCTCCACCAGTTTTGGACCGGATGAGCCCA  
GCTGGGTCCCTGCTGCCTACTTGAGAAAGTGGTACGCTGTACCCATACACCCGGCAGAAGGACAATGA  
GCTCTCTTTTCTGAAGGAACCGTCATCTGTGCTACTCGACGCTACTCAGATGGTGGTGTGAGGGTGTG  
AGCTCAGAGGGCACTGGATTCTTCCAGGGAATATGTGGAGCCAGCTG**TGA**

**ACGGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001013118  
**Insert Size:** 1104 bp



[View online »](#)

|                               |   |
|-------------------------------|---|
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).  |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <u><a href="#">NM_001013118.1</a></u> , <u><a href="#">NP_001013136.1</a></u>   |
| <b>RefSeq Size:</b>           | 1708 bp   |
| <b>RefSeq ORF:</b>            | 1104 bp   |
| <b>Locus ID:</b>              | 303476  |
| <b>Cytogenetics:</b>          | 10q26   |