

Product datasheet for **RN209435**

Pop5 (NM_001105752) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Pop5 (NM_001105752) Rat Untagged Clone
Tag: Tag Free
Symbol: Pop5
Synonyms: Rnasep3; rpP2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN209435 representing NM_001105752
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGTGCGGTTCAAGCACAGGTACCTTCTGTGCGAGCTGGTGTCCGAGGACCCGCGCTGCCGCTGAGCC
 TGGACGACCGGTGCTGGGCGGGCTCGTCCGGACACGATCGCCCGGTGCACGGGGCTTCGGCGCCG
 CGCTGCTCCGTGGGCTTGCAGTTAGATACCTCAATGCTTACACTGGAGTAGTGTCTTCGATGCCGA
 AAGGATTTCTACCAGCTCGTGTGGTCACTCTTCTTTCATCACGTGTTGGAGAACAAGGACACCGTT
 ACTCGTGTCTTCAACACATTACATGTGGGAGGTACCATTAGAAGCTGTCAGAAGTTCTGATTGAGTA
 CAACAGGAGACAGCTGTTGGTCTTGTGAGAAGTACACCGATGAAGGAGAGCGGGAAGCCATTAAGAAG
 TCTGTCTCAAGAAGCTGTCTGCTGGACAGAGAGTCGGTTGAAGAGCTTTCGGACAGTGTGGGAGGAGG
 CTGCAGAGGCAATGGAG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001105752
Insert Size: 510 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001105752.2, NP_001099222.2</u>
RefSeq Size:	1439 bp
RefSeq ORF:	510 bp
Locus ID:	117241
Cytogenetics:	12q16