

Product datasheet for **RN209453**

Olr60 (NM_001000748) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olr60 (NM_001000748) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Olr60
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN209453 representing NM_001000748 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGATCTCCTCACAGGCTTTCATCAACATTTCTTCTTCCAGCCACAGTCTTTCCTCATGATTGGCATCC
 CAGGACTAGAGGCTGTCCACGGCTGGATCTCCATCCCCTTCTCCTTCATGTACACTGTGGCACTCACC
 AAAGTGTCTCATCCTCCTGGCTGTGAGGAGGACCCATAGCCTACACCAGCCCATGTAATCTCTATCC
 ATGCTGGCCCTTAGTGATGTGGGCTCAGCTTGTCCACACTGCCCTCCACCCTGGCTGTGCTCTGGTTTG
 ACTATCGCTTCATCGGCTTCAATGCCTGTCTGATCCAGATGTTCTTCTACACTCCTTCTCCGTGTTGGA
 GTCCTCAGTGCTCTTGGCATGTCTTTGACCGATTTGTGGCTATCTCCAATCCACTGCGCTATGCATCT
 GTCCTCACTAACACTGTCTCATCATCAGGATCGGGGTGGCCATTGTAGCTAGGGCCACTCTGTCACTATTCC
 CAGTGCCCTTCTTGCTGAAACGACTGAAGTACTGCCCTGGCAAGATCCTCCTGTCTCACTCATTCTGTTT
 CCATGCTGATGTCTGAAGCTGGCTGTGCTGACATTACTGTCAACATCCTCTATGGACTCTATGTGGTT
 CTGTCCACAGTGGGTGTAGACTCCTCGCTTATTGTCTGCTACTCCTGATTCTTTACACAGTGATGG
 GGCTAGCCTCTCCAAGGAGCGCGTCCGACCCCTCAACACATGTGTTTCCCATATCTTGGCTGTCTGGT
 CTTTACATTCCAGTCATAGGTGTGTCCATGATCCACCGTTTGGAAAGCACCTACCTCACATTGTTTCAT
 GCCCTTGTGCGTACGTGTACCTTGTGGTACCTCCTGTGCTCAACCCCATCATCTATAGTGTCAAGTCCA
 AGCCCATCAGGGAGCCATGTTCTAGAGTTCTAAGAGGAAAAGACCAGACCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001000748
Insert Size:	963 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001000748.1, NP_001000748.1</u>
RefSeq Size:	963 bp
RefSeq ORF:	963 bp
Locus ID:	405017
Cytogenetics:	1q32
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]