

Product datasheet for RN209512

Npffr2 (NM_023980) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Npffr2 (NM_023980) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Npffr2
Synonyms:	Gpr74; Npff2; Npgpr
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN209512 representing NM_023980 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGCAAGAGATGGGACTCAAACCTCTTCAGGAAGCTGGGATCACATCTGGAGTGGCAATGACACACAGC
ATCCTTGGTATTCAGATATCAACATCACATACATGAACACTACTATCTCCACCAGCCCCACGTGACAGCTGT
CTTCATTAGCTCCTACTTCTGATCTTCTCTGTGCATGGTGGGAAACACTGTCGTTTGTCTTTGTGTA
ATAAGGAATAGGTACATGCACACGGTCACTAATTTCTTCATCTTCAACCTCGCAATAAGTGACTTACTGG
TTGGAATATTCTGCATGCCTATCACATTGCTGGACAACATCATAGCAGGATGGCCGTTTGAAGCAGCAT
GTGCAAGATCAGCGGGCTGGTGAAGGGATATCGGTTGCCGCTTCTGTCTTACCTTGGTTGCCATAGCC
GTAGACAGATTCCGGTGTGTGGTCTACCCCTTTAAGCCCAAGCTCACTGTCAAGACAGCCTTTGTCATGA
TCGTGATCATCTGGGCTGGCCATCACCATTATGACCCCATCTGCAATCATGTTACATGTACAGGAAGA
AAAATACTACCGTGTGAGGCTCAGCTCCCAATAAAAACCAGCACAGTCTACTGGTGTGGGAGGATTGG
CCAAACCAGGAAATGAGGAGGATCTACACCACCGTGTCTTTGCCACTATCTACCTGGCTCCACTCTCC
TCATTGTTATCATGTATGCAAGGATTGGGGCTTCCCTCTTCAAGACCTCAGCACACAGCAGGTAAGCA
GCGCCTGGAGCAGTGGCATGTATCCAAGAAGAAACAGAAGGTCAAGATGCTGCTGACTGTGGCCCTC
CTTTTCATCCTTTCTGGCTTCCCTGTGGACTCTGATGATGCTCAGACTATGCTGACCTGTACCTA
ACAAACTACGTGTCATCAATATTTATGTCTACCCTTTTGCCCACTGGCTCGCCTTCTGCAATAGCAGTGT
CAACCCCATCATTATGTTTTCTTTAATGAAAATTTTCGAGTGGTTTTCAAGATGCTTTCCAGTCTCTGC
CAAAAGAAAGTCAAACCCAGGAAGCCTATGGCCTAAGAGCTAAACGCAACCTGGACATAAACACATCTG
GCCTGTTGGTCCATGAACCTGCATCTCAAACCCAAGTGGGAAAACCTGGGATGTAGAAAAAGTGCAGA
CAATCCACACAGGAATCCTTGATGGAGGAAACGGGAGAAGCTACCAACAGTACTGAGACT**TAG**

ACGGTACGGCGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Chromatograms:	https://cdn.origene.com/chromatograms/ja2782_e01.zip
Restriction Sites:	Sgfl-Mlul
ACCN:	NM_023980
Insert Size:	1254 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_023980.1 , NP_076470.1
RefSeq Size:	1254 bp
RefSeq ORF:	1254 bp
Locus ID:	78964
UniProt ID:	Q9EQD2
Cytogenetics:	14p22
Gene Summary:	a receptor for neuropeptide FF (NPFF); may be involved in pain modulation and opiate tolerance [RGD, Feb 2006]