

Product datasheet for **RN209543**

Srxn1 (NM_001047858) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Tag: Tag Free

Symbol: Srxn1

Synonyms: Ab2-390; Npn3

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN209543 representing NM_001047858
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGTACAGCCCGGCCGGGTAAAGAGGCCGGGTGGCGCTGTCCGCCAAGGAAGCAATATGGGGCTACGCG
CGGGTGGAGCTCTCCGAGGGCGAGTGCAGGGCCCTGGGGCGCCGACGGTCAGGGACCCAGTGGCGCGCA
AGGAGGCAGCATCCACTCTGGCTGCATCGACACTGTGCACAACGTGCCAATCGCCGTGCTCATCCGGCCG
CTGCCGTGGTGTGGACCCGGTCAAGGTGCAGAGCCTCGTGGACACGATCCTGGAGGACCTGACAGTG
TGCCCCCATCGACGTCCTCTGGATCAAAGGGGCCCAAGGCGGTGACTACTACTCTTTGGGGGCTG
CCACCGCTATGCAGCCTACCAGCAGCTACAGAGAGAGACCATTCTGCCAAGCTCGTGAGGTCTACCCTT
TCAGACCTGAGGATGTACCTGGGTGCATCCACACCAGACTTGCA**GTAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001047858

Insert Size: 468 bp



OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001047858.3, NP_001041323.3</u>
RefSeq Size:	2878 bp
RefSeq ORF:	468 bp
Locus ID:	296271
Cytogenetics:	3q41