

Product datasheet for **RN209570**

Rwdd3 (NM_001128152) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rwdd3 (NM_001128152) Rat Untagged Clone
Tag: Tag Free
Symbol: Rwdd3
Synonyms: Rsume; X2cr1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN209570 representing NM_001128152
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGAGGAGATGAGGCAGAACTCTCGGCGCTGGCCGCGATCTTCTGCGGGCCCCACGAGTGGGAGA
TGTTGAGCTGTTAGAGACAGATGGGGCTGTCTTTAGAATTCACACAACAGCTGAAGGACTTGCGGGTGA
GGATGTGCCCTTAGAGTTGGCGTCCACTTACCAGCCGTTACCCATCGTGTCTACCCGGCATTTCGGTT
AACTCGGAACGTCTGACCAGGGCCAGTGTGTACCGTGAAGGAGAACTACTCGGAGAAGCCAGGAGGC
TTCTGTGAGAACCAGTGGTTCACGAGCTGGTCTTTGGACTCAGCAGAATCTCAGACACATCCTCAGCCA
GACAGAGACAGAGAGCAATGGGACGTGACTCTGCCAGAAAGCTCGACTGTGGACGGCGGACTGTGG
ATGACTCTTTTGCCTTAGATCAGATGAGAGCGAGGACTAAATATGTGAAAGTTGTGGAGAAGTGGGCTT
CAGAGTTAAGGCTGACAGGAAGACTGATGTTTATGGGTAAGTAAATGATACTGATTTTACTACAAGGAGACAG
AAGCAACATCAAGGAGTACCTGATTCTTCAGAAAACCTCAAAGTAGATGTGGACTCAAGTGGAAAGAAA
TGCAAAGAGAAAATGATTAGTGTACTGTCTGAAACAAAAGTACAAACGGAACACAAAAGGTTTCTGGCGT
TTGAAGTCAAAGAGTATTCAACTGGAGGAGTTACAAAAGGAATTTGAACTGCAGGACTTCAGGAGCT
CTTCTCAGAGTGTACTTGGCTGGTGAATGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001128152
Insert Size: 804 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001128152.1</u> , <u>NP_001121624.1</u>
RefSeq Size:	1044 bp
RefSeq ORF:	804 bp
Locus ID:	65026
UniProt ID:	<u>P0C7N0</u>
Cytogenetics:	2q41
Gene Summary:	expressed in GH3 pituitary cells [RGD, Feb 2006]