

Product datasheet for RN209702

Tsen34 (NM_001006968) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tsen34 (NM_001006968) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Tsen34
Synonyms:	LENG5; SEN34
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN209702 representing NM_001006968 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGGTGGTGGAGGTCGCTAATGGCCGCTCGCTGGTATGGGGAGCCGAGGCGGTGCAGGCGCTGCGGG
AGCGCCTGGGAGTCGGGGACGCACGGTGGGCGCCCTGCCCGCGGGCCCGCCAGAAGCTCGCGCTGGG
CCTTCCGCTTCTGCTGCTGCCTGAAGAAGCCCGCTCCTGGCGGAGATAGGCGCTGTGACCTTAGTCAGC
GCCCGCGTCCGGATCCCGCAACCATGGCTTGGCCCTAGCGTCGTTCAAACGCCAGCAAGAGCAGAGTT
TCCAGGATCAGAGCACTTGGCAGCAGAGGCCCGGAGACCCGGCGTCAGGAGCTTCTGGAGAAGATCGT
AGAAGGCCAGGCTGCCAAGAAGCAGAAGCTGGAACAGGATTCTGGGGCAGATGAAGAAGCCAAGAAGCC
GGTGAAGTGAAGCTACCAAGGAAGTGAGACCAGTGATGATGGTCAGGCTTCTGCGGAGCAGGAAGGAG
CAGATTCTTATCCCCCAACCAGGACCTTCAAATGGGGTGACTCCCTTGCCAGATCAGCCCTGCTTAT
CCAGCTGGCCACTGCCAGGCCTCGGCCTGTAAGGCTAAGCCTCTGGACTGGCGTGTGCAGTCCAAGAC
TGGCCCATGCTGGCCGCTCCTGCTCATGAGCTTCGCTACAGCATCTACAGAGACCTGTGGGAGAGAGGCT
TCTTCTCAGCGCAGCTGGGAAGTTTGGTGGTACTTCTGGTCTATCCTGGTATCCGCTGCGTTCCTCA
TGCTACTACATTGCTCAGTCTGGTCTGCTGAGGACCCATCCCACTTCAGGACCTGGTCTCTGCAGGC
CGCCTGGGAACCAAGTGTGAGAAAGACTTGCTTCTGCTCCCTCAGCCTGATGGGAAGGTGGTTTATA
CTTCCCTGCAGTGGGCCAGCCTGCAGTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001006968



[View online »](#)

Insert Size:	939 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001006968.1</u> , <u>NP_001006969.1</u>
RefSeq Size:	1369 bp
RefSeq ORF:	939 bp
Locus ID:	292534
Cytogenetics:	1q12