

Product datasheet for **RN209778**

Ccdc86 (NM_001006974) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ccdc86 (NM_001006974) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Ccdc86
Synonyms:	MGC94720
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN209778 representing NM_001006974 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGGTACGCCGCTGAGGCGCAGCCGGCGGCTGGAAGGCCTAAACCTCTGTCTCTGGAGAACCTTCCTG
ACCCCGAAGTCTCTCGGGCCAAACGGGCTCTTGTAGATTTCAAGTCGAACCCAGAAGAAACAAGAGAGCT
CGAGTCTCCGAGAGTTCCGCCACTGGACCTTGTGTCTCCCAACCTCAGCCGGAGACTAGCCCCGAATCA
CCCTGTCTAAACAGGATGCAGGCTTTGGGTCCCCACAAAGCAACCAGAGCCTCACCTGGGTCCCTTC
AACTTCATCAAGATCTGGGCCTGGATTACCTGCAGGACAGACAGAGTCGAATCCGGAATCCCCGCAACG
AGAGCAGTCATCCAAATTATCTCCACACAGGACTCTGAGGTAGCCACGCTAAGGAGGAGGTGATCCCA
AGATCTCCCGAGCCTTGCCCGGGTCAGCAAGCACCGGGTCCCGAGCCCTCCAGCCAGCACAGGAGCTGG
CATTCCAGGCACCTAGCTCCCCAGAGCGCCAGCTGGAGCCAAGCAAGCTACCTCCGGCTGGAGAGTCAGT
GACCGGAAGTCTCGACCTAAAGAAGCGAGTTATAGCTTCTCCCGAGGCCCGAGCTCCAAGAACTGAAG
GAGGAGCTTCTGTAAATCCGAAGGGGAAGCCCAATCCGGGCGTGTGTGGAAAGACCGCTCCAAGAAGA
GGTTCTCCAGATGGTTCAGGACAAGCCCTGCGTACCTCCTGGCAGAGGAAGATGAAGGAGCGACAGGA
GAGGAAGCTGGCTAAGGACTTCGCCCGGCACCTGGAGGAGGAGAAACAGAGGCGCCGGCAGGAGAAGAAA
GAACGCCGAGCTGAGAACCTCAGACGCCGTCTAGAAAACGAACGCAAGCTGAGATTGTCCAAGTGATCC
GGAACCTTGCCAAGCTCAAGAAAGCAAGAAGAAGCAGTTGCGCTCCATTCAAAGAGAGACACCTGGC
CTTGCTCCAGAAACAGCCACCCAGAGGCTGTAGCCAAGGTC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001006974
Insert Size:	1026 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001006974.2</u> , <u>NP_001006975.1</u>
RefSeq Size:	1026 bp
RefSeq ORF:	1026 bp
Locus ID:	293738
UniProt ID:	<u>Q5XIB5</u>
Cytogenetics:	1q43