

Product datasheet for **RN209833**

Rpl13a (NM_173340) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rpl13a (NM_173340) Rat Untagged Clone
Tag: Tag Free
Symbol: Rpl13a
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN209833 representing NM_173340
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGAGGGGAGGTTCTAGTATTGGATGGCCGGGCCATCTTCTGGGCCGCTGGCGGCCATTGTGG
CCAAGCAGGTACTGCTGGGCCGAAAGTGGTGGTTGTACGCTGTGAGGGCATCAACATTTCTGGAAATTT
CTACAGAAACAAGTTAAAGTATCTGGCCTTTCTCCGAAAGCGGATGAACACCAACCCGCTCGAGGCCCC
TACCACCTCCGAGCCCCAAGCCGATTTTTGGCGCACTGTGCGAGGCATTGCTCCGACAAAGACCAAAA
GAGGCCAGGCTGCCCTGGAACGCCTCAAGGTGTTGGATGGGATCCCTCCACCCTATGACAAGAAAAAGCG
GATGGTGGTCCCTGCTGCCCTCAAGGTTGTGCGGCTGAAGCCTACCAGAAAGTTTGCTTACCTGGGGCGT
CTGGCTCATGAGGTCGGGTGGAAGTACCAGGCAGTGACAGCTACTCTGGAGGAGAAACGGAAGGAAAAGG
CAAAGATCCATTACCGAAGAAGAAGCAGCTCTTGAGGCTAAGGAAACAGGCAGAAAAGAATGTGGAGAA
GAAAATCTGCAAGTTCACAGAGGTCTCAAGACCAATGGACTCTTGGT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_173340
Insert Size: 612 bp
OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_173340.2</u> , <u>NP_775462.2</u>
RefSeq Size:	690 bp
RefSeq ORF:	612 bp
Locus ID:	317646
Cytogenetics:	1q22
Gene Summary:	60S ribosomal subunit protein; mouse homolog is mutated in the tum-transplantation antigen P198 [RGD, Feb 2006]