

## Product datasheet for RN209878

### Tssk1b (NM\_001011900) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tssk1b (NM\_001011900) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Tssk1b  
**Synonyms:** Stk22a; Tssk1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN209878 representing NM\_001011900  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGC**C

ATGGATGACGCAGCCGCTCTCAAGCGACGAGGCTACATCATGGGGATAAACTTAGGAGAGGGCTCGTATG  
CCAAAGTCAAATCCGCTTACTCTGAGCGCTAAAGTTCAACGTGGCGGTCAAGATCATCGACCGCAAGAA  
AGCCCCACGGACTTCTGGAGAAATTCCTTCCCGGGAGATTGAGATTCTGGCCATGCTAAACCACCGC  
TCTATTGTCAAGACCTATGAGATCTTTGAGACGTCAAGTGGCAAGGTCTACATTGTGATGGAGCTTGGGG  
TCCAGGGTGACCTCCTTGAATTCATCAAAACCAGAGGAGCCCTGCAAGAGGATGATGCTCGAAAGAAGTT  
CCACCAGCTCTCCTCGGCATCAAGTACTGCCACGACCTGGATGTCGTCACCGAGACCTCAAGTGTGAG  
AACCTTCTGCTTGACAAGGACTTCAACATCAAAGTGTCTGACTTTGGCTTCTCCAAGCGCTGCTTGGGG  
ATGATAGTGGTCGGCTGACATTAAGCAAGACCTTCTGTGGGTGAGCAGCTTATGCGGCCCGGAGGTGCT  
GCAGGGCATCCCCACCAGCCAAAGGTGTATGACATCTGGAGCCTGGGTGTGATCCTCTACATCATGGTC  
TGTGGCTCCATGCCCTACGATGACTCCAACATCAAAAAGATGCTGCGCATCCAGAAGGAGCACCCGCTCA  
ACTTCCCACGCTCCAACACCTGACTGGTGAGTGCAAGGACCTCATCTACCGAATGTTGCAGCCAGATG  
CAACCGGAGGCTGCACATTGATGAGATCCTCAGCCACTGCTGGGTGCAGCCCAAGACTCGAGGCTGTCC  
TCTGGAGCCATCAACAAGGAAGGGGAAAGCTCCCGGGCCGCTGAGCCCTCATGGACCCCTGAACCTGGTG  
CTGATAAGAAGTCTGCCACCAAGTTGGAGCCTAAGGAAGAGGCACGAGCTGAGGCCCAATCTGAGACAAA  
ACCCAGGAAGATAATGTGACAGGTGGCCAGGCAGTCAGAGAACGTGGGCCTCAGCAGCGAACTGAACAGG  
GCCACAGAGGAAGGCATCCCCAACAGCCTTCAGAGACACATACT**TAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001011900



<b>Insert Size:</b>	1098 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001011900.1, NP_001011900.1</u>
<b>RefSeq Size:</b>	1447 bp
<b>RefSeq ORF:</b>	1098 bp
<b>Locus ID:</b>	288358
<b>Cytogenetics:</b>	11q23