

Product datasheet for **RN209897**

Spout1 (NM_001047110) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Spout1 (NM_001047110) Rat Untagged Clone
Tag: Tag Free
Symbol: Spout1
Synonyms: MGC112863
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN209897 representing NM_001047110
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTGTCCAGCTGGCCCCATCCTGCAGTACCTGGAGTGTCCACAGTAAGTGTGAGGTATCTGAGGAAAG
CGTTCTTCCCAAGCACCAGGATCTGCAGTTTGCAGGAATCCTCAACCCGTTGGACAGCCCTCACCACAT
GCGTCAGGATGAAGAATCTGAGTTCAGAGAAGGCGTCGTGGTGGACCGGCCACCAAGGCTGGCCACGGC
TCTTTGGTCAACTGTGGAATGAAGAAGGAAGTGAAGATTGATAAGAACTGGACCCTGGACTTCGGGTGA
CCGTGCGGCTGAACCAGCAGCAGCTCCCAGAATGCAAGACGTACAAAGGAACAGTCGTGTCGTACAGGA
CCCCCGCACAAAGCTGGTCTCTACTGGGGCTACTGTCCGGCTGGCCCTCTGCCTCAGTGTGTGTTT
GCCGAGGCTCCCTTCAGGATGGGTATGACCTGACCATTGGGACATCAGAACGTGGCTCCAGCGTGGCCCT
CTACCCAGCTGCCCAACTTCAGGCATGCTCTTGTGGTGTTCGGAGGCTCCAAGGCTAGAAGCTGGAGT
GGATGCTGACCCGAACCTTGAAGTGGCTGACCCAGTGTCTCTTCGACTTTTATGTCAACACATGCCCT
AGCCAGGGCAGTCGTACCATCCGCACCGAGGAGGCCATACTCATCTCACTGGCCGCTCTGCAGCCTGGCC
TCACCCAGGGAGGGCGCGACCCACTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001047110
Insert Size: 729 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001047110.1</u> , <u>NP_001040575.1</u>
RefSeq Size:	1623 bp
RefSeq ORF:	729 bp
Locus ID:	499770
Cytogenetics:	3p12