

## Product datasheet for **RN209935**

### **Ccdc94 (NM\_001109673) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ccdc94 (NM\_001109673) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Ccdc94  
**Synonyms:** MGC156684  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >RN209935 representing NM\_001109673  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTCGGAACGGAAAGTTTTAAACAAATACTACCCACCTGACTTTGACCCATCGAAGATCCCGAACTCA  
AGTTACCCAAGGACCGGCAGTACGTGGTGCGGCTGATGGCCCCATTCAACATGAGGTGTAAGACATGTGG  
AGAATACATCTACAAGGGGAAGAAGTTCAATGCACGAAAGGAAACAGTCCAGAATGAGGCATACCTGGC  
TTGCCAATCTTCCGTTTCTATATCAAGTGCACAGCTGCTTGCCGAGATCACCTTTAAGACAGACCCAG  
AGAACACAGACTACACCATGGAGCATGGGGCCACGCGCAACTCCAGGCTGAGAAGCTTCTGGAGGAGGA  
GGAAAAGCGCATGCAGAAGGAGCGTGAGGATGAGGAGCTGAACAACCCCATGAAGGTCCTAGAGAACCGT  
ACAAAGGACTCTAAGCTGGAGATGGAAGTGCTAGAGAACCTACAAGAGCTAAAAGACCTGAACCAGCGGC  
AGGCCACGTGGACTTTGAGGCCATGCTGCTACAACACCGCCTATCTCAGGAGCAGTGGCAACAGCAGCA  
GGAAGAGGAAGACGAGCGTGAGACAGTGGCCTTGCTGGAGGAGGCTCGTTATCGAAGGCTTCTGGATGAC  
TCTGACTCAGAGGATGAGGCTCCACCTTCCAACCACAGGCAGCTGCAAGACCCAACCCACAGCCATTC  
TTGATGAGGTACCAAAGACCAAGAGGAAGGCAGAGGGCCTGTGCAGGAAAGCCAGTTGGCAGGCTTGGT  
TATACCAAAGAAGGCAAAGACGGAAGCCAACCGGTTTCAGAGCAGGTCGGGTGCCACCACCTGGTGCC  
CCCGAGGCAGGAAGGTGGCCAACCCCTGCACCCAGACACCTGGTGCCTCCTCCCTGAGCCAGCTCGGTG  
CATATGGGGACAGTGAGGACAGTGACAGC**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001109673



<b>Insert Size:</b>	942 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001109673.1</a></u> , <u><a href="#">NP_001103143.1</a></u>
<b>RefSeq Size:</b>	2096 bp
<b>RefSeq ORF:</b>	942 bp
<b>Locus ID:</b>	501285
<b>Cytogenetics:</b>	9q12