

Product datasheet for **RN210013**

Serhl2 (NM_001130579) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Serhl2 (NM_001130579) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Serhl2
Synonyms:	RGD1563386
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN210013 representing NM_001130579 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGGTTTGCTCTCAGAGTTGAAGTTGGCTGTGCCCTGGGGCCACATTGCTGTCAAAGTCTGGGGCTTGC
AGAAGAACCCTCCAGTTCTGTGCCTGCATGGCTGGTTGGACAATGCCAACTCCTTTGACAAGCTCATTCC
TTTTCTTCCAAAAGACTTTTGTATGTGGCCATGGATTCGGAGGCCATGGGCTTTCAACCCACTACAGC
CCAGGTCTCCCTACTACCACCATAACTTCGTGAGTGAGGTCCGAAGGGTGGTGTGAGCTTTTAAGTGGA
CTCGGCTCTCCCTCCTGGGCCACAGTTTGGTGGCGTTGTGGGTGGCTTGTTCGCCTGCATGTTCCCTGA
GATGGTGGACAACTGATCTTGCTGGACTCGACGCCACTTCTCATGGACTTGAATGAAGTGAGAAACATA
ATGACCTACAGGCGGAAGAACATAGAGCAGACACTGAAGGTGGAAGACTCCCAGAAGCCCCGGGCATCT
TCAGCCCAGAGGAGATGCTGCATGAGCTCCTGACTAAGAACAGCCACTTGAATGAAGACTGTGGGAACT
CCTCCTGCAGAGGGGAACCACGAAGGTGGCTGAAGGTCTGGTGTGAACAGGGACCAGAGGCTCTTTGG
CCACAGTACAGCTTTGACTTCATGGGCAAGGAGCTGACTATGCACTCTTTGAGGAGGCTGCAGGCCAGCG
TGTTAATAATCAAAGCACTGGACGGGTACTACGACGTGAGACGGGAGAACGATTTGAATAAGCGTCCTT
CCTCCTCATGCTGGACATTCTACGATCCACCCTAAAGGAGCGGTTCCAGTTTGTGGAAATCCCAGGCAAT
CACTACATCCACATGAACAAGCCCCAGATCGTGGCTGACATCATCCGCTCCTTCTACAGGGCCACCGA
GGACAACCCCTTCCAGGCTG**AG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001130579


[View online »](#)

Insert Size:	933 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001130579.1</u> , <u>NP_001124051.1</u>
RefSeq Size:	2493 bp
RefSeq ORF:	933 bp
Locus ID:	500911
Cytogenetics:	7q34