

Product datasheet for **RN210113**

Psme2 (NM_017257) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Psme2 (NM_017257) Rat Untagged Clone
Tag: Tag Free
Symbol: Psme2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN210113 representing NM_017257
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCAAGCCTTGTGGGTTTCGCTTGAGTGGAGAAGCCCGTAAACAGGTGGATGCCTTCAGACAAAATC
TTTTCCAGGAGGCTGAGGACTTCCTTTGCACTTCTTGCCACGGAAAATCATATCCCTGAGTCAGCTCTT
GCAAGAGGATTCCTCAATGTGGCTGACCTCTCCTCCCTCCGGGCTCCTCTGGACATCCCTATCCCAGAT
CCCCACCCAAGGATGATGAGATGGAAACCGAGCAGGAGAAGAAAGAAGTCCCTAAGTGTGGCTTCTCC
CTGGGAATGAGAAGCTTCTGGCCCTGCTTGCTTGGTTAAGCCAGAAGTCTGGACTCTCAAAGAGAAGTG
CATTCTGGTAATCACATGGATCCAGCACCTGATCCCCAAGATTGAGGACGGAAATGATTTTGGGGTGGCA
ATTCAGGAGAAGGTGCTGGAGAGAGTGAATGCTGTCAAGACCAAAGTGGAAAGCCTTCCAGACAGCCATCT
CCAAGTATTTCTCAGAGCGTGGGATGCTGTAGCCAAGGCCTCCAAGGACACCCATGTAATGGATTACCG
GGCCTTGGTGCACGAGCGAGATGAAGCAGCCTATGGGGCGCTTAGGGCCATGGTGTGGACCTGAGAGCC
TTCTATGCTGAGCTTCATCACATCATCAGCAGCAACCTAGAGAAAATTGTCAACCCAAAGGGTGAAGAAA
AGCCATCTATGTACT**GTA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul
ACCN: NM_017257
Insert Size: 717 bp



| | |
|-------------------------------|---|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_017257.2 , NP_058953.1 |
| RefSeq Size: | 1366 bp |
| RefSeq ORF: | 717 bp |
| Locus ID: | 29614 |
| UniProt ID: | Q63798 |
| Cytogenetics: | 15p13 |
| Gene Summary: | may play a role in proteasome activation [RGD, Feb 2006] |