

Product datasheet for **RN210207**

Alg13 (NM_001013951) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Alg13 (NM_001013951) Rat Untagged Clone
Tag: Tag Free
Symbol: Alg13
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN210207 representing NM_001013951
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGAGAGCGTTTGTGACCGTCGGGACCACCAGTTTCGACGACCTCGTCGCTCGTGTGGTCGCCAACG
ATACGGTTCAGATACTCAAGAGTCTGGTTACAACCACCTTGCTCTACAAATTGGTAGAGGCACTGTGGT
GCCTGAGCCATTCACTAGTCTGAGCCATTCACTCTGGATGTTTACAGGTATAAGGAGTCTTTGAAAGAAGAC
CTTCAGCAAGCTGATCTTGTCATTAGCCACGCAGGTGCAGGAAGCTGTTGGAGAGTCTAGAAAAAGGCA
AGCCACTTGTGGTAGTTGTAATGAAAAGTTAATGAACAATCATCAATTTGAACTGGCAAAGCAGTTACA
CAAAGAAGGCCATCTCTTCTATTGTACTTGACGATGCTTCTGGGCTGTACAGTCAATGGATTTATCA
ACACTGAAATGTTATCTCCTGGCCAGCCAGAAAAATTTCTGCATTTTTGGATAAAGTTGTTGGATTAC
AAAA**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001013951
Insert Size: 498 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001013951.1](#), [NP_001013973.1](#)

RefSeq Size: 1625 bp

RefSeq ORF: 498 bp

Locus ID: 300284

UniProt ID: [Q5I0K7](#)

Cytogenetics: Xq34

Gene Summary: May be involved in protein N-glycosylation, second step of the dolichol-linked oligosaccharide pathway.[UniProtKB/Swiss-Prot Function]