

Product datasheet for **RN210386**

Dusp22 (NM_001108412) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dusp22 (NM_001108412) Rat Untagged Clone
Tag: Tag Free
Symbol: Dusp22
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN210386 representing NM_001108412
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGGAGTGGGATGAGCCAGATCCTGCCGGCCTGTACATTGGCAACTCAAAGATGCAAGAGATGCAG
AACAAATTGAGCAGGAACAAGGTGACACACATTCTTTCTGTGCACGACACTGCCAGGCCATGCTGGAGGG
AGTTAAATACCTGTGTATTCCAGCAGCAGACTCACCGTCTCAAACTGACAAGACATTTCAAAGAAAGT
ATTAATTCATTGATGAGTGCCGACTCCAGGGTGAGGGCTGTCTGTACATTGCCTGGCTGGGGTCTCCA
GGAGTGTGACACTGGTGATCGCATAACATCATGACTGTCACCGACTTTGGCTGGGAAGAGGCCCTTACACAC
CGTTCTGTGACGGGAGGTCTGTGCCAACCCCACTGGGCTTCAAAGGCAGCTGCAGGAGTTTGAGAAA
CACGAAGTGCGGCAGTATCGGCAATGGCTGAGAGAAGAATATGGAGAGAACCCTTGCGGGATGCAGAAG
AAGCCAAAAGCATTCTGGCTGCCCGGGAATTCTGAAGTACTGGGCTTTCTCAGAAGACTG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001108412
Insert Size: 555 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001108412.1, NP_001101882.1

RefSeq Size: 1158 bp

RefSeq ORF: 555 bp

Locus ID: 361242

Cytogenetics: 17p12