

Product datasheet for **RN210388**

Meis1 (NM_001134702) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Meis1 (NM_001134702) Rat Untagged Clone
Tag: Tag Free
Symbol: Meis1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN210388 representing NM_001134702
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGCAAAGGTACGACGATCTACCCATTATGGGGGCATGGATGGAGTAGGCATCCCCTCCACGATG
ATGGGGACCCGCATGCAGCCAGGTCTATGCAGCCGGTCCACCACCTGAACCACGGGCCTCCTCTGCACTC
TCATCAGTACCCGCACACAGCTCATACCAACGCCATGGCCCCAGCATGGGCTCCTCAGTCAATGACGCT
TTAAAGAGAGATAAAGATGCCATTTATGGACACCCCTCTCCCTCTCTTAGCACTGATTTTTGAGAAAT
GTGAATTAGCTACTTGTACCCCGCGAGCCGGGGTGGCGGGCGGGACGTCTGCTCGTCAGAGTCATT
CAATGAAGATATAGCGGTGTTGCGCAAACAGATTCGTGCAGAAAACTCTATTCTTTCTAATCCAGAA
CTGGATAACTTGATGATCCAAGCCATACAAGTATTAAGGTTTCATCTGTTGGAATTAGAGAAGGTACACG
AATTATGTGACAATTTCTGCCACCGGTATATTAGCTGTTGAAAGGGAAAATGCCTATCGATTTGGTGAT
AGATGATAGAGAAGGCGGATCAAATCAGACAGTGAAGATGTAACAAGAGCAGCAAATCTAACTGACCAG
CCCTCTTGGAATAGAGACCATGATGACACGGCATCCACTCGTTCAGGAGGAACCCCGGGCCCTTCCAGCG
GTGGCCCACTTCACACAGTGGGATAACAGCAGTGAAGGTGATGGCTTGGACAACAGTGTAGCTTC
CCCCAGCAGGTCAGATGATGACCCTGATAAGGACAAAAGCGTCACAAAAGCGTGGCATCTTTCCC
AAAGTAGCCACCAATATCATGAGGGCGTGGCTGTTCCAGCATCTAACACCCCTTACCCTTCTGAAGAGC
AGAAAAAGCAGTTGGCACAAGATACGGGACTCACCATCCTTCAAGTGAACAATTGGTTTATTAATGCGCG
GAGAAGAATAGTGAGCCATGATAGATCAGTCCAACCGAGCAGTCAAGGACACCTTATAACCCC
GATGGACAGCCAAATGGGAGGTTTGTAAATGGACGGTCAGCAGCACATGGGCATCAGAGCGCCAGGACCTA
TGAGTGAATGGGCATGAATATGGGCATGGAGGGCAGTGGCACTACATG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001134702



[View online »](#)

Insert Size:	1173 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001134702.1, NP_001128174.1</u>
RefSeq Size:	2506 bp
RefSeq ORF:	1173 bp
Locus ID:	686117
Cytogenetics:	14q22